VDW.GOLD® RECIPROC®

de	Gebrauchsanweisung	2
en	Operating Manual	37
fr	Mode d'emploi	71
it	Manuale d'uso	106
es	Manual de instrucciones	141

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**Congratulations on your purchase of the
VDW.GOLD® RECIPROC® endo motor.**

Endodontic products and services are the sole focus of VDW GmbH's business operations. Constantly remaining in close cooperation with international universities and endodontists enables us to develop innovative concepts, products and systems which make dentists' work simpler, easier and more effective.

By having purchased the VDW.GOLD® RECIPROC® endo motor, you have acquired a product which has been developed and tested with the utmost care, meeting even the highest expectations in terms of functionality and operation.

VDW GmbH reserves the right to edit the information and data in this manual without prior notice.

The manual is available in additional languages upon request.

This manual has been created with the greatest care, but even with all our efforts, errors can never be entirely excluded. Suggestions for improvements are welcome anytime. Please contact VDW GmbH with any related comments or suggestions you may have.

VDW GmbH
Bayerwaldstr. 15
81737 Munich
Germany

Phone +49 89 62734-0
Fax +49 89 62734-304
info@vdw-dental.com
www.vdw-dental.com

Table of Contents


1. Symbols Used	39	7.5 Selecting File Systems/Files	53
1.1 In these Instructions for Use	39	7.5.1 Changing Torque and Rotational Speed (only in rotary mode)	53
1.2 On Packaging, Device and Components	39	7.5.2 Dr's Choice (only in rotary mode)	54
		7.5.3 Rotary Mode	54
2. Indications for Use	40	7.5.4 Reciprocation Mode	54
		7.5.5 Integrated Apex Locator	55
3. Contraindications	40	7.5.6 The Apex Locator LED Display	55
		7.5.7 Combined Length Determination	56
4. Warnings	40	7.5.8 Separate Length Determination (Working Length Determination with Hand Instrument)	57
		7.5.9 Valuable Information and Suggestions on Precise Length Determination	57
5. Precautions	42	7.6 Default Parameters	58
6. Adverse Reactions	42	8. Maintenance, Cleaning, Disinfection and Sterilisation	58
7. Step-by-Step Instructions	42	8.1 Maintenance	58
7.1 Standard Components	43	8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664)	60
7.2 Getting Started	44	8.2.1 Pre-Treatment	60
7.2.1 Preparation	44	8.2.2 Manual Cleaning and Disinfection	61
7.2.2 Overview of Connections	44	8.2.3 Inspection/Maintenance	61
7.2.3 Power Supply	45	8.2.4 Packing	61
7.2.4 Micromotor and Foot Pedal	45	8.2.5 Sterilisation	62
7.2.5 Apex Locator	45	8.2.6 Storage	62
7.2.6 Battery	46	8.2.7 Material Durability	62
7.3 User Interface	47	9. Technical Data	63
7.3.1 Keypad	47	10. Troubleshooting	64
7.3.2 Display	48	11. Warranty	69
7.3.3 Foot Pedal	48	12. Dr's Choice Tables	70
7.3.4 Audible Signals	49		
7.3.5 File Library	50		
7.4 Operation	50		
7.4.1 Power On, Standby Mode, Power Off	50		
7.4.2 CAL: Calibration	51		
7.4.3 ASR: Automatic Stop Reverse	51		
7.4.4 ANA: Anatomy Function	52		
7.4.5 MENU: Device Settings	52		
		Appendix	
		Electromagnetic Emissions and Immunity (English)	177

1. Symbols Used

1.1 In these Instructions for Use

	<p>Failure to follow all instructions correctly during use can damage the product or endanger users/patients.</p>		<p>Additional information; explanation of operation and performance.</p>
---	---	---	--

1.2 On Packaging, Device and Components

	<p>Serial number</p>		<p>Consult instructions for use</p>
	<p>Manufacturer</p>		<p>Not to be disposed of with ordinary waste (in accordance with Directive 2002/96/EEC regarding the separate collection of electrical and electronic devices).</p>
	<p>Date of manufacture</p>		<p>Foot pedal connection</p>
	<p>Protection Class II Device</p>		<p>Lip clip cable connection</p>
	<p>Type BF applied part</p>		<p>File clip cable connection</p>
	<p>Caution</p>		<p>Direct current (power-supply connection)</p>
	<p>CE marking</p>		<p>Autoclavable (note the temperature indicated on the symbol)</p>
	<p>Temperature limitation</p>		<p>GOST symbol: product meets Russian security standards (GOST-R).</p>
	<p>Humidity limitation</p>		<p>Catalogue number (reorder number)</p>
	<p>Keep dry!</p>		<p>Only use the device with the power supply provided.</p>
	<p>Fragile!</p>		

2. Indications for Use

FOR DENTAL USE ONLY!

VDW.GOLD® RECIPROC® is a medical device in accordance with the Medical Device Directive 93/42/EEC, amended with the Directive 2007/47/EC. The endo motor is exclusively designed for use by dentists with dental root canal instruments in continuous rotation and in reciprocating movement with integrated apex locator (endometric length determination).

This device may only be used by qualified dental professionals in hospital settings, clinics or dental practices.

3. Contraindications

Do not use the device for implants or other non-endodontic dental procedures.

4. Warnings

This chapter contains a description of severe adverse reactions and potential safety hazards for the product or the user/patient. Please read the following warnings before use.

WARNINGS

Users

• VDW.GOLD® RECIPROC® is designed for use in dental care and may only be operated by trained, qualified professionals, e.g. dentists.

Ambient Conditions

- VDW.GOLD® RECIPROC® should not be placed in damp areas or anywhere where the device may come into contact with any type of liquids.
- Do not expose the device to any direct or indirect sources of heat. The device must be operated and stored in a safe environment.
- Do not use the device in the presence of free oxygen, anaesthetic substances or flammable products. The device must be operated and stored in a safe environment.

- The device requires special safety precautions with regard to electromagnetic compatibility (EMC), and the EMC information contained in this manual must be strictly adhered to during installation and operation. It is especially important not to use the device near fluorescent lamps, radio transmitters or remote controls, portable or mobile radio frequency communications equipment.

- In order to avoid possible danger from electromagnetic interference, do not operate any other electromedical or other electrical devices near the VDW.GOLD® RECIPROC®. The electromagnetic radiation emitted by this device is below the recommended limits set by the applicable regulations (DIN EN 60601-1-2:2007).

- Do not pinch the cables on the micro-motor, the lip clip, the file clip, the foot pedal or the charger, and do not pull directly on the cables when plugging in or unplugging components. Make sure cables are no obstacle and do not limit mobility.

- The use of accessories other than those specified, with the exception of components sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the device.

Power Supply

- Only use the battery charger provided.
- The battery charger must be connected to a power supply measuring 100-240 V (+/- 10%), 47-63 Hz.
- Turn the VDW.GOLD® RECIPROC® off before disconnecting the charger.
- Do not position the battery charger in a way that it is difficult to disconnect the device from the mains.

Battery

- Always work on battery power to ensure long battery life. It is recommended that you only recharge the battery after it has completely discharged.
- If the battery LED light begins to flash red during use, connect the device to the battery charger immediately in order to prevent it from shutting down.

- Never open the device to change the battery, as this could cause a short-circuit. Opening the device voids the warranty. Batteries may only be replaced at your service centre.

- Should any liquid exit from the VDW.GOLD® RECIPROC® that might be attributed to a leaking battery, discontinue use of the device immediately and send it to your service centre for battery replacement.

During Use

- Dental dams and gloves are mandatory when using the VDW.GOLD® RECIPROC®.

- Weak auxiliary currents flow through the electrodes during endometric length determination. Maximum values for the VDW.GOLD® RECIPROC® are well under the limits set by IEC 60601-1:2005. Nonetheless, in rare cases, these currents can cause patients to feel pain. Should this occur, discontinue the length determination process.

- Should any irregularities occur when using the device, switch the micromotor off and contact your service centre.

Components and Accessories

- Only use the included lip clip cable with the ferrite ring.

- Only use the VDW 6:1 contra-angle with the VDW.GOLD® RECIPROC®. The precision of the torque, the rotational speed and the length determination are only ensured if the VDW 6:1 contra-angle is used.

The micromotor changes rotational speed from the minimum to the maximum value during calibration. Do not insert any files during calibration.

- Always calibrate after the VDW 6:1 contra-angle has been cleaned and sterilized, or at least once per week (see separate operating instructions for the VDW 6:1 contra-angle).

- Precision of the micromotor's motion is only guaranteed if the original VDW 6:1 contra-angle is used and properly maintained. When maintaining the contra-angle, please prevent lubricant from penetrating the micromotor. Lubricant can damage the micromotor and negatively affect its ability to operate safely. Before maintenance, read the separate operating instructions for the VDW 6:1 contra-angle as well as chapter

8.1 Maintenance in this manual.

- Do not insert any foreign objects into the micromotor shaft.

- Excessive application of force can cause the micromotor to overheat. If the micromotor overheats too often or fails to cool down after overheating, contact your service centre.

- Before starting the micromotor, make sure the settings are correct. The file system indicated on the display must always correspond to the file being used. This is extremely important in order to prevent reciprocating files from being used in rotary mode and vice versa.

- Only use endodontic NiTi instruments with metal shafts for combined length determination.

- Follow the instructions given in the file manufacturer's manual when using endodontic NiTi instruments. These are not included in delivery.

- Never use rotating files in reciprocating mode. Check the mode shown on the display before use.

- Never use reciprocating files in rotary mode. Check the mode shown on the display before use.

- Torque and rotational-speed values can be changed by the file manufacturers without prior notice. The preinstalled values must therefore be checked in the library before use. The values shown on the display are only precise and reliable in combination with a correctly maintained and lubricated VDW 6:1 contra-angle.

- VDW.GOLD® RECIPROC® may not be connected to or used in combination with other devices or systems. It may not be used as a component for another device or system. Never connect external PC storage media (e.g. hard drives) to the VDW.GOLD® RECIPROC® via the USB port. The USB port is exclusively intended for use by authorized personnel for purposes of maintenance or software updates. The manufacturer assumes no responsibility for any accidents, damage to the device, bodily injury or other problems resulting from non-compliance with this limitation.

Care and Transportation

• The VDW.GOLD® RECIPROC®'s components are not delivered disinfected or sterilised: components such as the central unit, the micromotor and the micromotor cable as well as the lip clip and file clip cables must be disinfected before the first use and after every additional use. The VDW 6:1 contra-angle, the lip clip and the file clip (without cables) must be sterilised before the first use and after every additional use!

• Never put the micromotor or other accessory parts in an autoclave or an ultrasonic bath. Do not sterilise any VDW.GOLD® RECIPROC® components (with the exception of the lip clip, the file clip (not including cables) and the VDW 6:1 contra-angle). See separate operating instructions on how to sterilise the VDW 6:1 contra-angle.

• The plastic enclosure is not sealed. Do not use any liquids or sprays directly on the console, especially on the display or near the electrical sockets.

Repair and Malfunctions

• Do not use the VDW.GOLD® RECIPROC® if you suspect it is damaged or defective.

• Repairing, changing or modifying the VDW.GOLD® RECIPROC® without the manufacturer's prior approval is not permitted. VDW GmbH assumes no responsibility if it has been changed or modified. Should a malfunction occur, contact your service centre instead of having repairs done by an unauthorized person.

Disposal

• Disposing of the VDW.GOLD® RECIPROC® in the normal household waste is not permitted. Dispose of the device in accordance with Directive 2002/96/EEC regarding the separate collection of electric/electronic devices. For further details, please contact VDW GmbH directly.

• For information on the proper disposal of the endodontic NiTi instruments, please refer to the manufacturer's instructions. These are not included in delivery.

5. Precautions

Read all safety information carefully before use. These precautionary measures help you use the product safely and thus prevent injury to both you and others. Keep this user's manual in a safe place for future reference. This manual must be included whenever the system is sold or passed on so that the new owner can observe these warnings and take the appropriate precautions.

All special precautions to be taken before, during and after operation can be found in Chapter 4 Warnings. The manufacturer assumes no responsibility in case of:

• Use of the VDW.GOLD® RECIPROC® for purposes other than those specifically described in this manual.

• Performance of work or repairs by anyone not authorized by the manufacturer or by the importer in the manufacturer's name.

• Use of non-original components or components other than those described in Chapter 7.1 Standard Components.

• File breakage in the root canal during clinical use of the VDW.GOLD® RECIPROC®.

• Connection of the micromotor to a power source which does not conform to IEC 364 norms.

• Damage to device or accessories resulting from sterilisation: VDW.GOLD® RECIPROC® components cannot be sterilised (with the exception of the lip clip, the file clip (not including cables) and the VDW 6:1 contra-angle. See also the separate operating instructions for information on sterilising the VDW 6:1 contra-angle.)

Should you have any questions, please contact your supplier or the VDW GmbH service centre.

6. Adverse Reactions

There are no known adverse reactions.

7. Step-by-Step Instructions

This chapter contains all necessary information regarding the commissioning and operation of the VDW.GOLD® RECIPROC®.

NOTE

Before using the device for the first time, familiarize yourself with Chapter 4 Warnings which contains all special precautions to be taken before, during and after operation.

7.1 Standard Components

The VDW.GOLD® RECIPROC® is shipped with the components listed below:


Fig. 1 VDW.GOLD® RECIPROC®

1) Central unit
2) Micromotor with cable and connector (cable length 1.8 m)
3) Reduction gear contra-angle: VDW 6:1 contra-angle (packed in a separate shipping box with accompanying operating instructions)
4) 2 protective silicon covers for the VDW 6:1 contra-angle
5) 2 lip clips
6) Lip clip cable (with ferrite ring); cable length 1.7 m
7) 2 file clips
8) File clip cable (cable length 1.7 m)
9) Foot pedal with cable (length 1.7 m)
10) Battery charger (with ferrite ring) with interchangeable EU, UK, USA, and AUS plug adaptors (cable length 1.8 m)
VDW.GOLD® RECIPROC® user manual

The VDW 6:1 contra-angle is packed separately and delivered in a shipping box along with the VDW.GOLD® RECIPROC®. Note the separate operating instructions for the VDW 6.1 contra-angle.

WARNING

The standard components are delivered neither sterilised nor disinfected! Kindly note the respective information in Chapter 8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664) before operating the device.

7.2 Getting Started

7.2.1 Preparation

1. Carefully remove the central unit and the accessories from the packaging and set them on a flat surface in the desired location. Keep the environmental specifications listed below in mind when choosing this location.

2. After opening the package and before beginning installation, check the device for possible damage and for missing parts. Notify your supplier within 24 hours after receiving the device about any missing parts or damage in transit.

3. Check that the device serial number on the bottom of the central unit matches the one indicated on the package.

4. Check that the serial number of the VDW 6:1 contra-angle matches the one indicated on its box.

5. Sterilise the following components in the autoclave:

- VDW 6:1 contra-angle (consult the separate operating instructions)
- Lip clip and file clip (not including cables), see Chapter 8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664).

6. Disinfect the front of the central unit, the micromotor, the micromotor cable, and the lip clip and file clip cables; see Chapter 8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664).

WARNING

Should any liquid exit the device, discontinue use immediately and send the device to your service centre.

Environmental Conditions Required for Operation

Use: in closed rooms

Ambient temperature: 15 °C - 42 °C

Relative humidity: < 80 %

The original packaging materials can be stored and shipped at temperatures between -20 °C and +50 °C with relative humidity levels of 20 - 90 %.

WARNING

Do not set the device in a damp place or anywhere where it will come into constant contact with fluids of any kind.

7.2.2 Overview of Connections


Fig. 2 Front Connection


Fig. 3 Rear Connections

Front View

A) Connection socket for the micromotor

Rear View

B) Connection socket for the charger

C) Connection socket for the lip clip cable

D) Connection socket for the file clip cable

E) Connection socket for the foot pedal

F) USB port for updates (only for service personnel)

NOTE

The USB port (Fig. 3, F) is used to install software updates or perform maintenance on the device. This port is only to be used by authorized service personnel and is covered by a label reading "ONLY FOR SERVICE DO NOT REMOVE" on delivery. The label may only be removed by service personnel.

7.2.3 Power Supply

1. Choose the plug adaptor appropriate to your electricity supply (see Fig. 4).
2. Place the appropriate adaptor onto the two contacts (Fig. 4, G) on the charger and push it towards the lock button until it clicks into place. Push the lock button to change the plug adaptor.


Fig. 4 Charger and plug adaptor

3. Connect the charger cable into connection socket (Fig. 3, B) at the back of the device.
4. Use the power cord to connect the device to the mains. The battery will begin charging.

NOTE

Charge the battery completely before using the device for the first time. For information on proper battery care, see Chapter 7.2.6 Battery.

7.2.4 Micromotor and Foot Pedal

5. Insert the plug of the micromotor cable into connection socket (Fig. 2, A) at the front of the central unit.

NOTE

The connection socket (Fig. 2, A) is a push-pull connector. Align the red dot on the plug vertically and position it so that it can be inserted into the socket. Do not screw the plug into the socket. Pull on the metal end of the cable to unplug it. Do not twist the plug when removing.

6. Remove the black spray adaptor from the VDW 6:1 contra-angle.
7. Insert the VDW 6:1 contra-angle into the micromotor. (See also the separate operating instructions for the VDW 6:1 contra-angle).
8. Insert the foot pedal cable into connection socket (Fig. 3, E) at the back of the central unit.

7.2.5 Apex Locator

9. Insert the lip clip into the holder of the lip clip cable (with ferrite ring).
10. Insert the plug of the lip clip cable into connection socket (Fig. 3, C) at the back of the central unit.
11. Connect the file clip to the file clip cable.
12. Insert the plug of the file clip cable into connection socket (Fig. 3, D) at the back of the central unit.

7.2.6 Battery

The VDW.GOLD® RECIPROC® runs on a nickel-metal hydride (NiMH) battery. This chapter provides all important information regarding proper handling of the battery.

Battery LED Indicator

The battery LED indicator uses three different colours to display current battery status:

Green: Indicates a battery capacity between 20 and 100 %.

Flashing red: The battery needs to be recharged. Battery power will be exhausted within the next few minutes.

Before the motor shuts off, an audible alarm will sound for 22 seconds at two different frequencies and the following message will appear:

BATTERY

Flashing orange: The charger has been connected correctly and the battery is being charged. The following message will appear on the display:

VDW.GOLD RECIPROC
BATTERY CHARGING

VDW.GOLD® RECIPROC® can be used normally while charging without significantly lengthening the time needed to recharge. The charger is strong enough to supply the micromotor with electricity directly.

Once the battery is fully recharged, the battery LED light will turn green again.

! WARNINGS

- If the battery LED light begins to flash red during operation, connect the device to the charger immediately, as the micromotor could otherwise shut down during the procedure.
- The battery warning is a HIGH priority alarm signal.

Operation and Recharge Times

Recharging the battery after it has been discharged completely can take up to 3 hours.

If the battery has been correctly charged, it should have an operation time of at least 2 hours before needing to be charged again.

Prolonging Battery Life

In order to ensure long battery life, it is recommended to always work on battery power and to only recharge once the battery has discharged completely.

For optimal battery performance, the battery should be replaced at your service centre once every 3 years.

! WARNINGS

- Never open the device yourself in order to change the battery, as this can cause it to short-circuit. Opening the device voids the warranty. The battery should thus only be replaced by personnel at your service centre.
- Should the device begin to leak fluid which might have something to do with a leaking battery, discontinue use immediately and send the device to your service centre for battery replacement.

7.3 User Interface

7.3.1 Keypad


Fig. 5 Keypad

1) ON/OFF

Turns the device on and off.

2) Battery LED indicator

Indicates the current battery status (see details under 7.2.6 Battery).

3) Apex Status LED indicator

Indicates the current status of the apex locator:

Green LED: Apex locator on

Red LED: Apex locator off

4) Apex Locator LED indicator

Reference scale for root canal length determination (see details on this in Chapter 7.5.6 The Apex Locator LED Display).

5+6) |◀ / ▶|

You can use the right/left keys in the bottom row of the display to:

- page through the File System (**sys**), File (**file**), Torque (**gcm**) and Rotational Speed (**rpm**) fields. The active field is indicated by an arrow (only available in rotary mode).

- select individual menu (**MENU**) functions. Current settings are indicated by arrows.

7+8) +/-

You can use the +/- keys to:

- scroll up/down through file systems or files and change individual torque or rotational speed values in all rotating file systems.

- Select settings on individual menu **MENU** functions. Current settings are indicated by arrows.

9) ✓ Confirm

- Saves changes to torque or rotational-speed settings in all modifiable systems.

- Also restores standard parameters when pressed simultaneously with the **ON/OFF** key (see details in Chapter 7.6 Default Parameters).

- The **✓ Confirm** key can also be used to start the micromotor without using the foot pedal.

10) ANA

Activates (green LED on) or deactivates (LED off) the ANA function (see details in Chapter 7.4.4 ANA: Anatomy Function).

11) MENU

You can change the following settings using the **MENU** key (see details in Chapter 7.4.5 MENU: Device Settings)

APEX STOP

LANG (LANGUAGE)

SYSTEM

SOUND

RECIPROC Reverse Function (RECIPROC-REVERSE): This function is activated at the time of delivery. It can be deactivated (OFF) or reactivated (ON) (see details in Chapter 7.5.4 Reciprocation Mode).

12) ASR

You can use the ASR key to select various ASR functions (details in Chapter 7.4.3 ASR: Automatic Stop Reverse):

Green LED: Auto-stop rotational direction reversal in rotary mode once preset torque value has been reached.

Red LED: Reversal of rotational direction without torque limits in rotary mode – controlled by the foot pedal.

LED off: In rotary mode, the micromotor stops once the preset torque value has been reached.

ASR is deactivated in reciprocating mode (LED off).

13) CAL

Calibrates the VDW 6:1 contra-angle in order to ensure the precision of its torque values following replacement or lubrication (see details in Chapter 7.4.2 CAL: Calibration).

7.3.2 Display

When the device is switched on, all LEDs light up and a welcome text is displayed. During subsequent uses, the display shows the last first used before the device was powered off.

When you bring the device out of standby mode (by pressing the foot pedal or any key on the keypad), the display will also show the file last used before standby mode was activated.

Top Line:

Indicates the selected mode – either rotary or reciprocation – by displaying either ROTARY or RECIPROCATION.

Bottom Line in Rotary Mode:

Displays the File System (**sys**), File (**file**), Torque (**gcm**) and Rotational Speed (**rpm**) fields:

sys

Displays the selected rotary file system (e.g. MTWO for Mtwo®, FM for FlexMaster®, DR's for DR'S CHOICE, etc.)

file

Displays the selected file.

gcm

Shows the set torque limit value (deactivated in reciprocation mode). Torque is indicated in g/cm (grams of force per centimetre: 1 gcm = 0.0981 Nmm).

If the standard torque settings have been changed, a **!** symbol will be displayed in front of the new value.

rpm

Shows the rotational speed of the instrument (deactivated in reciprocation mode). Rotational speed is indicated in rpm (revolutions per minute).

If the standard rotational speed settings have been changed, a **!** will be displayed in front of the new value.

Bottom Line in Reciprocation Mode:

Displays the reciprocation system in use:

Unlike in rotary mode, individual files cannot be selected in reciprocation mode, as all motor settings apply to the entire system, indicated by RECIPROC ALL.

7.3.3 Foot Pedal

The micromotor can be started in two different ways:

- The foot pedal is pressed and the micromotor runs for as long as it remains pressed, or
- The micromotor is started by pressing the **✓ Confirm** key for 1.5 seconds. The micromotor can be stopped by pressing the foot pedal or any key on the keypad. If the micromotor is started using the **✓ Confirm** key, it will automatically shut off after 5 minutes of inactivity in order to save electricity.

7.3.4 Audible Signals

All audible signals are activated on delivery:

Components	Audible Signal	Setting
Keypad		
Keys	Confirmation sound when keys on the keypad are pressed.	Settings cannot be changed
Motor		
Rotary mode	Warning signal when torque exceeds approximately 75% of the preset value in rotary mode.	Volume control 0-3 Can be set under MENU --> SOUND --> MOTOR
	Intermittent signal during continuous counter-clockwise rotation.	Volume control 0-3 Can be set under MENU --> SOUND --> MOTOR
Reciprocation mode	Warning signal when the entire length of the instrument is affected by an increased stress. See Chapter 7.5.4 Reciprocation Mode for details.	Volume control 0-3 Can be set under MENU --> SOUND --> MOTOR
Apex Locator (accompanying warning signal for the apex locator LED display)		
3 blue LEDs	Sound sequence (Median coronal canal section)	Volume control 0-3 Can be set under MENU --> SOUND --> APEX
3 green LEDs/ 1 orange LED	Increasing signal frequency (sound sequence). When the orange LED has been reached, a constant beep will sound (apical constriction area to the apical foramen)	Volume control 0-3 Can be set under MENU --> SOUND --> APEX
1 red LED	Very loud warning signal (file surpassed the apical foramen (over-instrumentation))	Volume control 0-3 Can be set under MENU --> SOUND --> APEX
Battery		
Battery status	Loud warning signal which sounds before the device shuts down due to low battery power. The battery warning is a HIGH priority alarm signal.	Volume control 0-3 Can be set under MENU --> SOUND --> MOTOR

7.3.5 File Library

The device contains a file library with the following preset NiTi systems:

Reciprocating Systems

- RECIPROC®
- WaveOne™

Rotary Systems

- Mtwo® (MTWO)
- FlexMaster® (FM)
- DR'S CHOICE (DR'S)
- ProFile® (PF)
- System GT® (GT)
- ProTaper® Universal (PTU)
- Hero® (HERO)
- K3™ (K3)
- Race™ (RACE)
- FlexMaster® Retreatment (FMR)
- Lentulo (LENTU)
- Gates Glidden (GATE)

WARNINGS

- *The file system indicated on the display must always correspond to the file being used. This is extremely important in order to avoid using reciprocating files in rotary mode or vice versa.*
- *The values displayed are only precise and reliable in combination with a correctly-maintained VDW 6:1 contra-angle.*
- *Follow the file manufacturer's instructions when using the endodontic files.*
- *The manufacturer reserves the right to update the file library and the systems contained therein. The preinstalled values must therefore be checked in the library before use.*
- *Never use rotating files in reciprocating mode. Check the mode shown on the display before use.*

- *Never use reciprocating files in rotary mode. Check the mode shown on the display before use.*

7.4 Operation

7.4.1 Power On, Standby Mode, Power Off

Power On

Press the ON/OFF key. The LEDs will light up briefly as their functionality is tested. A welcome screen will display the current software version:

```
VDW.GOLD RECIPROC
SOFTWARE X.X
```

In subsequent sessions, the display will show the last file in use before the device was powered off.

When you bring the device out of standby mode, the display will show the file last used before standby mode was activated.

```
ROTARY
> MTWO 10/04 120 280
   sys   file   gcm   rpm
```

Standby Mode

In order to save battery power, the device will automatically go into standby mode after 10 minutes of inactivity. The display will shut off.

To leave standby mode, press the foot pedal or any key on the keypad. The device will switch on and return to the function last in use before it went into standby mode.

After 30 minutes in standby mode, the device will shut off completely in order to save battery power. It can be switched on again using the ON/OFF key.

Power Off

Press the ON/OFF key to shut off the device. The battery LED light will no longer be illuminated unless it is currently being recharged.

7.4.2 CAL: Calibration

Calibrating the micromotor automatically sets the rotational speed in order to ensure torque precision. Calibration is necessary:

- after replacing the micromotor,
- after replacing, cleaning or autoclaving the VDW 6:1 contra-angle, but at least once per week (see also the separate operating instructions for the VDW 6:1 contra-angle).

Starting Calibration

1. Put the VDW 6:1 contra-angle on the micromotor (without file).
2. Press the **CAL** key. The calibration process will begin, and the following message will be displayed.

The micromotor will automatically measure the inertia of the contra-angle by starting at minimum rotational speed and accelerating to maximum speed.

3. The procedure stops automatically as soon as calibration is complete (all black bars are visible on the display).


WARNINGS

- *Never press the CAL key while the device is in use.*
- *Calibrate the device without files. The micromotor accelerates from minimum to maximum rotational speed during calibration.*
- *Use only the VDW 6:1 contra-angle with the VDW.GOLD® RECIPROC®. The precision of the torque, the rotational speed and the length determination are only ensured if the VDW 6:1 contra-angle is used.*


Aborting Calibration

The calibration procedure can be halted at any time by pressing the foot pedal or any key on the keypad. The display will then read:

Possible Error Messages

If calibration cannot be correctly performed or completed, the display will read:

Details on possible error messages such as ERROR 1 or ERROR 2 can be found in Chapter 10: TROUBLESHOOTING.

7.4.3 ASR: Automatic Stop Reverse

NOTE

The ASR function (automatic stop and reversal of rotational direction) can only be activated in rotary mode. In reciprocation mode, the ASR function is auto-matically deactivated (LED off).

ASR is active (green LED) in rotary mode on delivery.

Selection of Different ASR Functions in Rotary Mode

You can select different ASR functions by pressing the ASR key. The LED light will change colour depending on the function selected:

Green: When the preset torque value has been reached, the micromotor will automatically rotate in the opposite direction (counter-clockwise). As soon as the file stops meeting resistance, the micromotor will automatically begin running forwards (clockwise) again.

Red: The micromotor will automatically rotate in the opposite direction (counter-clockwise) without torque control.

OFF: When the preset torque value has been reached, the micromotor will stop automatically. Pressing the foot pedal again will start the motor in the opposite direction (counter-clockwise). As soon as the file stops meeting resistance, the micromotor will automatically begin running forwards (clockwise) again.

NOTE

ASR functions can be set for each file in rotary mode. As soon as you select a different file, the standard settings (green LED) will be restored.

Audible Warning Signals sound:

- whenever torque exceeds approximately 75% of the preset value;
- when the micromotor is running in reverse.

7.4.4 ANA: Anatomy Function

NOTES

- *The ANA function can only be activated in rotary mode.*
- *In reciprocation mode, the ANA function is automatically deactivated (LED off).*

The anatomy function (ANA) automatically reduces the selected instrument's torque in the VDW 6:1 contra-angle. This safety feature is generally recommended for new users and in cases of complicated canal anatomy.

This function is inactive on delivery.

Activating and Deactivating the ANA Function

To activate the ANA function, press the ANA key (green LED light). To deactivate the ANA function, press the ANA key again (LED off).

7.4.5 MENU: Device Settings

You can change the following settings in the menu:

APEX STOP: Activate (ON) or deactivate (OFF) apical auto-stop (see details in Chapter 7.5.7 Combined Length).

LANGUAGE (LANG): Select the menu language (GERMAN/ENGLISH).

SYSTEM: On device delivery, the major reciprocating and rotary file systems available in the VDW.GOLD® RECIPROC® file library are shown (see details in Chapter 7.3.5 File Library). File systems can be shown or hidden via the menu according to personal reference. You can hide (OFF) any file systems you do not need under the SYSTEM menu item. These file systems will then no longer be displayed. To show a file system again, open the SYSTEM menu item and set the desired file system back to ON.

SOUND: Set the volume of the warning signals for the apex locator LED display (0-3) and the motor (0-3) (see details in Chapter 7.3.4 Audible Signals).

RECIPROC REVERSE Function:

This function is activated on delivery. It can be deactivated (OFF) and reactivated (ON) as needed (see details in chapter 7.5.4 Reciprocation Mode).

Changing Settings

1. Press the MENU key.
2. Use the +/- keys to select the appropriate menu item (scroll up/down).
3. Use the |◀/▶| keys to select the desired setting, and then change it using the +/- keys.
4. Save your new settings by pressing ✓ Confirm.

7.5 Selecting File Systems/Files

The display's standard setting is to show the last file system or corresponding last file used. To make a new selection:

1. Use the **|◀/▶|** keys to select the **sys** field (arrow):

2. Use the **+/-** keys to select the desired file system (scroll up/down).

3. Use the **|◀/▶|** keys to select the file field (arrow):

4. Use the **+/-** keys to select the desired file (scroll up/down).

5. It is not necessary to confirm the setting. The desired file or file system is now active.

7.5.1 Changing Torque and Rotational Speed (only in rotary mode)

In rotary mode, the preset torque and rotational speed values for all file systems can be changed individually.

NOTE

This only applies in rotary mode. Torque and rotational speed cannot be changed in reciprocation mode.

Changing Torque:

1. Select the desired rotary file as described above (arrow indicates active selection).
2. Use the **|◀/▶|** keys to select the **gcm** field (arrow):

3. Use the **+/-** keys to select the desired torque (scroll up/down). As soon as you have changed the torque value, it will flash. The available torque settings range from 20 to 500 gcm in steps of 10 gcm.

4. Save your new setting by pressing the **✓ Confirm** key once. The torque field will stop flashing and will be marked with a **⌘** symbol.

If the setting is not saved by pressing the **✓ Confirm** key, it cannot be used. The micromotor will not start.

Changing Rotational Speed:

1. Select the desired rotary file as described above (arrow indicates active selection).

2. Use the **|◀/▶|** keys to select the **rpm** field (arrow):

3. Use the **+/-** keys to select the desired rotational speed (scroll up/down). As soon as you have changed the rotational speed value, it will flash. The available rotational speed values range from 200 to 500 rpm in steps of 10 rpm, and from 500 to 2000 rpm in steps of 100 rpm.

The exception to this is Dr's Choice, for which you can set the rotational speed as high as 3200 rpm (up to 500 rpm in steps of 10, starting at 500 rpm in steps of 100 rpm).

4. Save your new setting by pressing the **✓ Confirm** key once. The rotational speed field will stop flashing and will be marked with a **⌘** symbol. If the setting is not saved by pressing the **✓ Confirm** key, it cannot be used. The micromotor will not start.

7.5.2 Dr's Choice (only in rotary mode)

You can use the Dr's Choice Individual Program to create your own file sequence. This enables you to manage your own series of instruments independently of file manufacturers or of recommended sequences.

The device is delivered with 15 standard torque and rotational speed values. You can find details as well as the "Factory Default Values in Dr's Choice Program" in Chapter 12 Dr's Choice Tables.

To change these settings individually, simply "overwrite" them as described in Chapter 7.5.1 Changing Torque and Rotational Speed (only in rotary mode). Use the table in Chapter 12 Dr's Choice Tables to record your custom settings.

To restore the standard settings, see Chapter 7.6 Default Parameters.

7.5.3 Rotary Mode

When a file system has been selected in rotary mode, the first line of the display will read:

ROTARY			
>	MTWO 10/04	120	280
	sys	file	gcm rpm

The second line will automatically display the first file in the sequence.

In rotary mode, torque and rotational speed values can be changed individually for all file systems. See Chapter 7.5.1 Changing Torque and Rotational Speed (only in rotary mode) for a description of this process.

WARNING

Never use reciprocating files in rotary mode. Check the mode indicated on the display before use!

NOTE

The ASR function can only be activated in rotary mode. The ASR function is active on delivery (green LED). Press the ASR key to select different functions (see details in Chapter 7.4.3 ASR: Automatic Stop Reverse).

7.5.4 Reciprocation Mode

RECIPROC® and WaveOne™ instruments are specially designed for use in reciprocation, whereby the instrument is driven first in a cutting direction and then reverses to release the instrument. The angles of reciprocation are precise and specific to the design of the instrument.

WARNING

Never use rotating files in reciprocating mode. Check the mode indicated on the display before use!

NOTES

- The RECIPROC® and WaveOne™ reciprocating files can be distinguished from rotary files by their special form: The shaft has a coloured plastic ring and their cutting edges are inverted.
- The ASR function, which is designed for constantly rotating files, is deactivated in reciprocation mode.

Working in Reciprocation Mode

When one of the reciprocating files mentioned above is selected, the first line of the display will read:

RECIPROCATION	
>	RECIPROC ALL

Unlike in rotary mode, individual files cannot be selected in reciprocation mode, as all motor settings apply to the entire system, indicated by RECIPROC ALL.

RECIPROC Reverse Function (for VDW RECIPROC® Instruments)

This factory preset function (see details in Chapter 7.4.5 MENU: Device Settings) enables more comfortable work and uses audible signals to alert the user to the following events or necessary actions:

- A low, fast beeping sound indicates increased stress along the entire length of the instrument. It is therefore necessary to enlarge the canal with lateral brushing motions. This will reduce the stress on the instrument.
- A higher, slower beeping sound indicates an even higher stress on the instrument. The micro-motor will therefore switch automatically to continuous clockwise rotation in order to release the instrument. Once the foot pedal has been released and pressed again, the micromotor will change back into reciprocating mode. It is now essential to brush laterally, to enlarge the root canal.

NOTES

Use an ISO size 10 C-PILOT® File to check that the canal is not blocked. Before brushing, repeat the following steps:

- Clean the instrument in the Interim Stand.
- Irrigate the canal.
- Use an ISO size 10 C-PILOT® File to check that the canal is not blocked.

If the RECIPROC Reverse Function has been deactivated (see details in Chapter 7.4.5: MENU: DEVICE SETTINGS), the motor will stop as soon as the factory-preset maximum torque has been reached. If this happens, remove the file from the root canal, clean the instrument and start over. However, we recommend that you always work with the RECIPROC REVERSE function, as it is more comfortable.

7.5.5 Integrated Apex Locator

The VDW.GOLD® RECIPROC® features an integrated apex locator for length determination of the root canal.

The apex locator can be used in 2 ways:

Combined Length Determination: The working length is determined while the root canal is being prepared. The micromotor and the apex locator are thus active simultaneously (use of the VDW 6:1 contra-angle and the lip clip).

Separate Length Determination: The working length is determined manually (without the micromotor) using the file clip and the lip clip.

7.5.6 The Apex Locator LED Display

As you are determining the length of the root canal, the apex locator LED display shows you the current position of the tip of the instrument.


Fig. 6 Apex Locator LED Display

The following areas are identified using a colour scale (LEDs):

- 1) **LED display:** 3 blue LEDs
Area: Median coronal canal section
Warning Signal: Slow sound sequence
- 2/3) **LED display:** 3 green LEDs/ 1 orange LED
Area: Apical section: Apical constriction area to apical foramen
Warning Signal: Increasing signal frequency (sound sequence). A constant beep will begin sounding when the orange LED has been reached.
- 4) **LED display:** 1 red LED
Area: Apical foramen surpassed (overinstrumentation)
Warning Signal: Very loud warning signal

NOTES

- The working length is set at the 3rd green LED.
- A single LED does NOT represent the working length in mm.

WARNING

Weak auxiliary currents flow through the electrodes during endometric length determination. Maximum values for the VDW.GOLD® RECIPROC® are well below the limits set by IEC 60601-1:2005. Nonetheless, in rare cases, these currents can cause patients to feel pain. Should this occur, discontinue the length determination.

7.5.7 Combined Length Determination

For a combined length determination:

1. Put the protective silicon cover onto the VDW 6:1 contra-angle.
2. Select a file (see details in Chapter 7.5 Selecting File Systems/Files).
3. To turn on the apex locator, touch lip clip with the file for about 3 seconds. **Do not touch the stopper!** See Fig.7.


Fig. 7 Activating the Apex Locator with Lip Clip and File on the VDW 6.1 contra-angle

When the Apex Status LED indicator lights up green, combined length determination has been activated, and the display will read:

APEX
LOCATOR ON

4. Hang the lip clip inside the patient's cheek (we recommend using the cheek opposite the tooth being treated).
5. Start the micromotor. The apex locator LED indicator shows you the current position of the tip of the instrument (see Chapter 7.5.6 The Apex Locator LED Display for details).
6. To interrupt or deactivate combined length determination, remove the lip clip from the patient's mouth or take the instrument out of the tooth.

Apical Auto-Stop

During combined length determination, apical auto-stop can be switched on and off as needed in both rotary and reciprocation modes.

If the function is on, auto-stop will automatically halt the instrument when the apex has been reached. Releasing and re-pressing the foot pedal will automatically turn the file in the other direction, then switch to the original direction.

Apical auto-stop is deactivated on delivery.

To activate/deactivate apical auto-stop:

1. Press the MENU key.
2. Use the +/- keys to select the menu item APEX STOP (scroll up/down).
3. Use the ►| key to select the ON/OFF field and change it to the desired setting (deactivate OFF/ activate ON) using the +/- keys.

MENU
APEX STOP >ON

4. Save your new settings by pressing ✓ Confirm.

7.5.8 Separate Length Determination (Working Length Determination with Hand Instrument)

You can use the file clip to determine the working length manually (without using the micromotor). Proceed as follows:

1. Select the desired hand instrument. Keep in mind: For precise results when doing separate length determination, the files used should be of a size appropriate to that of the apex. Using files that are too small will interfere with length determination.
2. To turn the apex locator on, attach the file clip for about 3 seconds to the lip clip. See Fig. 8.


Fig. 8 Activating the Apex Locator with Lip Clip and File clip

3. Insert the hand instrument into the separate file clip.
4. Hang the lip clip inside the patient's cheek (we recommend using the cheek opposite the tooth being treated).

When the Apex Status LED indicator lights up green, separate length determination has been activated, and the display will read:

**APEX
LOCATOR ON**

5. Bring the hand instrument into the root canal. The apex locator LED indicator shows you the current position of the tip of the instrument (see Chapter 7.5.6 The Apex Locator LED Display for details).
6. To interrupt or deactivate separate length determination, remove the lip clip from the patient's mouth or take the instrument out of the tooth.

WARNING

Weak patient auxiliary currents flow through the electrodes during endometric length determination. Maximum values for the VDW.GOLD® RECIPROC® are well below the limits set by IEC 60601-1:2005. Nonetheless, in rare cases, these currents can cause patients to feel pain. Should this occur, discontinue the length determination.

7.5.9 Valuable Information and Suggestions on Precise Length Determination

- Precise combined length determination is only guaranteed when the original VDW 6:1 contra-angle is used.
- Use only endodontic NiTi files with metal shafts for combined length determination.
- Use gloves and a dental dam to isolate the tooth.
- Dry the access cavity using an air blower or a cotton ball.
- Avoid direct contact between the file in the contra-angle and the oral mucosa, as stray current can disrupt the determination process.
- Use the protective silicon cover on the VDW 6:1 contra-angle.


NOTE

- *If error messages are displayed (e.g. APEX LOCATOR ERROR 3) or other problems occur, see Chapter 10 Troubleshooting.*

Endometric Length Determination and X-Ray Technology

Because X-ray measurement technology depicts the 3-dimensional root canal system in only 2 dimensions, the measurement results it provides may differ from those obtained endometrically. This is not an indication that the VDW.GOLD® RECIPROC® is unreliable, nor does it mean that the X-rays are inaccurate.

The differing lengths are caused by differences among canal anatomies. The actual apical foramen can be localised elsewhere from the radiological apex of the root.

If the canal is curved, the X-ray may show a shorter working length than the result determined using the VDW.GOLD® RECIPROC®.

7.6 Default Parameters

To restore the original standard parameters:

- Make sure the charger is not connected.
- Turn the device off.
- Press the **✓ Confirm** and **ON/OFF** keys at the same time. The device will power on and the display will read:

**DEFAULT PARAMETERS
LOADING**

NOTES

- *Keep in mind that all individual settings, including those for the Dr's Choice Program, will be erased if you restore the factory default settings.*
- *Settings regarding the audible signals are set to "2" after loading the default settings.*

8. Maintenance, Cleaning, Disinfection and Sterilisation

Maintaining proper hygiene when working with the VDW.GOLD® RECIPROC® is only possible with regular device care. Follow the maintenance, cleaning, disinfection and sterilisation instructions in the following chapters to ensure this.

8.1 Maintenance

NOTE

Service and repair work should only be done by factory-trained service personnel.

Components	Routine Maintenance
Cable	<p>Inspect the various cables of charger, micromotor, lip clip, file clip and foot pedal at least once every six months.</p> <p>Contact your service centre if any of the cable coatings show signs of wear and tear.</p>
VDW 6:1 contra-angle	<p>The contra-angle needs to be lubricated after cleaning and disinfection, but before sterilisation. Follow the instructions given in the separate operating instructions of the VDW 6.1 contra-angle, as well as the following notes:</p> <ul style="list-style-type: none"> • When lubricating the contra-angle, make sure no lubricant penetrates the micromotor. • If the contra-angle is being lubricated manually, use pressurized air to remove excess lubricant (blow for about 5 sec.) before re-applying the contra-angle back onto the micromotor. Calibrate the contra-angle after lubrication. • If the contra-angle is being lubricated automatically using a maintenance or lubrication device, follow the device manufacturer's instructions carefully, and make sure no excess lubricant remains in the contra-angle. • Do not lubricate the micromotor under any circumstances. Lubricant can contaminate the micromotor and create a safety hazard. Doing this invalidates the warranty.
Central unit	<p>Check whether the central unit is leaking fluids or emitting smoke. If so, unplug the device immediately from the power source and contact your service centre.</p>
Battery	<p>For optimal performance, the battery should be replaced at your service centre once every 3 years.</p> <p>Never open the device yourself in order to change the battery, as this can cause a short-circuit. Opening the device voids the warranty. The battery should thus only be replaced by personnel at your service centre.</p>

8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664)

Clean the cables and the outer surfaces of the device using a paper towel or soft cotton cloth lightly moistened with an aldehyde-free cleaning and disinfectant (bactericidal and fungicidal) solution, e.g. Mikrozid AF Liquid or Minuten Spray Classic.

Remove the lip clip and the file clip from the cables before cleaning. The lip clip and file clip need to be cleaned, disinfected and sterilised before each use. This includes first use. Thorough cleaning and disinfection are crucial for effective sterilisation. Follow the special instructions given in Chapter 8.2 Cleaning, Disinfection and Sterilisation (in Accordance with DIN EN ISO 17664). Consult the user's manuals of the devices used in your practice as well.

Ensure that only validated cleaning/disinfection and sterilisation methods are used, that the devices (disinfector, steriliser) are maintained and inspected at regular intervals, and that the validated parameters are preserved throughout every cycle.

Always observe the applicable government hygiene standards and your practice's or clinic's hygiene regulations as well. This applies especially to guidelines for effective prion inactivation.

For your own safety, always wear gloves when handling contaminated accessory parts.

WARNINGS

- *The cables cannot be autoclaved.*
- *Using substances other than those described above can cause damage to the device and the accessory parts.*
- *Do not use hot-air or radiation sterilisation, and do not sterilise using formaldehyde, ethylene oxide or plasma.*
- *The plastic casing is not sealed. Do not use any liquids or sprays directly on the console, especially on the display or near the electrical sockets.*

8.2.1 Pre-Treatment

1. Pulp and dentine residues must be removed from the accessory parts immediately (within 2 hours or less). Do not allow them to dry! After using the accessory parts on patients, put the parts directly into a bowl filled with an appropriate cleaning and disinfectant solution for cleaning, preliminary disinfection and temporary storage (for max. 2 hours).

2. Afterwards, clean the accessory parts under running water or in a disinfectant solution to remove all visible contaminants. The disinfectant should be aldehyde-free (aldehyde sets blood stains), tested for effectiveness (e.g. VAH/DGHM- or FDA-approved or CE certified), appropriate for disinfecting the accessories and compatible with the accessory parts (see Chapter 8.2.7 Material Durability).

Use only clean, soft brushes or cloths when removing contaminants manually. Do not use metal brushes or steel wool.

When cleaning the file clip, open and close it five times in order to clean the inner parts more effectively. Please note that disinfectants used during pre-treatment are only for personal protection and are not a substitute for post-cleaning disinfection. However, pre-treatment is still important.

WARNING

Do not use any automatic processes or ultrasonic baths to clean or disinfect the accessories.

8.2.2 Manual Cleaning and Disinfection

When selecting cleaning and disinfectant solutions, make sure that:

- they are intended for use in cleaning or disinfecting instruments,
- the disinfectant's effectiveness has been tested (e.g. VAH/DGHM- or FDA-approved or CE-certified) and that it is compatible with the cleaning solution,
- the chemicals used are compatible with the accessory parts (see Chapter 8.2.7 Material Durability).

Combined cleaning/disinfectant solutions should only be used if the instruments are only slightly contaminated (no visible contamination).

Follow the instructions provided by the manufacturers of the cleaning and disinfectant solutions regarding concentration, application time and rinsing.

Only use freshly-prepared solutions made with water which is either sterile or contains only low levels of bacteria (< 10 cfu/m) and endotoxins (< 0.25 EU/ml, e.g. purified water (PW/HPW)), as well as filtered, oil-free air for drying.

Step-by-Step Procedure

Cleaning

1. Place the pre-cleaned accessory parts in the cleaning bath for the stipulated application time; ensure that they are covered sufficiently (if necessary, use a soft brush to carefully brush them off). When cleaning the file clip, press it closed and release it five times in order to clean the inner parts more effectively.

2. Then take the instruments out of the cleaning bath and rinse them thoroughly with water at least three times for one minute each; press and release the file clip five times when rinsing.

Disinfection

1. After cleaning and inspecting the accessory parts, place them in the disinfection bath for the application time stipulated; make sure that the accessory parts are sufficiently covered. When disinfecting the file clip, open and close it five times in order to disinfect the inner parts more effectively.

2. Then take the instruments out of the disinfection bath and rinse them thoroughly with water at least five times for one minute each; open and close the file clip five times when rinsing.

3. Inspect, dry and pack the accessory parts as quickly as possible after taking them out of the disinfectant (see Chapters 8.2.3 Inspection/Maintenance and 8.2.4 Packing). Make sure that the accessory parts do not come into direct contact with one another!

8.2.3 Inspection/Maintenance

Check all accessory parts after cleaning or cleaning/disinfection. Dispose of defective accessory parts immediately. Defects may include:

- Warped plastic
- Corrosion

Accessory parts which are still contaminated must be cleaned and disinfected again. Maintenance is not required. Do not use instrument oil!

8.2.4 Packing

Pack the accessory parts in disposable sterilisation packaging (individual single-use packaging) which meets the following requirements:

- In compliance with DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilisation (temperature-safe to at least 137 °C (279 °F), sufficiently steam-permeable)

8.2.5 Sterilisation

The sterilisation procedure only applies to the lip clip and the file clip.

WARNING

Never put the micromotor or other accessory parts in an autoclave or an ultrasonic bath. Do not sterilise any VDW.GOLD® RECIPROC® components (with the exception of the lip clip, the file clip (not including cables) and the VDW 6:1 contra-angle. See also the separate operating instructions on sterilising the VDW 6:1 contra-angle.

Only use the sterilisation methods given below; other sterilisation procedures are not permissible.

- Steam sterilisation
- Fractional vacuum/pre-vacuum procedures (at least three vacuum cycles) or gravity sterilisation (product must be sufficiently dry). Gravity procedures are less effective and should only be used if fractional vacuum procedures are unavailable.
- Steam sterilisers in compliance with DIN EN 13060 or DIN EN 285
- Sterilisation validation must be done in accordance with DIN EN ISO 17665 (valid installation and operation qualifications (IQ and OQ) as well as product specific performance qualifications (PQ)).
- Maximum sterilisation temperature 134 °C (273 °F); plus tolerance as specified in ISO DIN EN ISO 17665.
- Sterilisation time (application time at sterilisation temperature) at least 18 min. at 134 °C (273 °F).

WARNINGS

- *Rapid-sterilisation processes and sterilisation processes with unpackaged accessory parts are not permitted.*
- *Do not use hot-air sterilisation, radiation sterilisation or sterilisation involving formaldehyde, ethylene oxide or plasma.*

8.2.6 Storage

After sterilisation, the instruments must be stored dry and dust-free in the sterilisation packaging.

8.2.7 Material Durability

When selecting cleaning and disinfection products, make sure that you do not use phenol, strong acids, strong aldehyde disinfectants or anti-corrosive solutions.

The material is temperature-resistant up to 137 °C / 279 °F (maximum application temperature).

9. Technical Data

MANUFACTURER	VDW GmbH, Bayerwaldstr. 15 81737 Munich, Germany
MODEL	VDW.GOLD® RECIPROC®
MEASUREMENTS	209 × 89 × 93 mm
MATERIALS	Console housing: PC/ABS Micromotor: Aluminium
WEIGHT	1.1 kg
POWER SUPPLY	NiMH Battery, 2000 mAh, 6 V
MAINS POWER SUPPLY	100-240 V
VOLTAGE FLUCTUATION	max. ± 10%
FREQUENCY	47-63 Hz
MAINS POWER SUPPLY – CURRENT CONSUMPTION NOMINAL VALUE	2.5 A
TORQUE RANGE	20-500 gcm (≈ 0.2-5.0 Ncm); +/- 30 % +/- 10 gcm in rotary mode
ROTATIONAL SPEED IN MICROMOTOR SHAFT	1200-19200 rpm (+/- 20 %) in rotary mode
ELECTRICAL PROTECTION CLASS	CLASS II
APPLIED PART	BF (contra-angle, lip clip, file clip)
SECURITY LEVEL NEAR FLAMMABLE ANAESTHETIC GAS MIXTURES OR OXYGEN	Not suitable for use near flammable anaesthetic gas mixtures or oxygen
OPERATIONAL MODE	Rotary and Reciprocation Modes
ENVIRONMENTAL CONDITIONS REQUIRED FOR OPERATION	+15 °C /+42 °C; RH: < 80%
MEDICAL PRODUCT CLASSIFICATION	Class IIa, Appendix IX, Rule IX, 93/42/EEC
CENTRAL UNIT AND MICROMOTOR	IP20
FOOT PEDAL	IPX1
TRANSPORT AND STORAGE CONDITIONS	-20 °C /+50 °C; RH: 20-90 %

10. Troubleshooting

Should the VDW.GOLD® RECIPROC® not function properly, it is not necessary indicative of a device malfunction. Use the information in the following table to check the device and rule out operating errors or anatomical or other irregularities.

Should the information in the following table not solve the problem, contact your service centre or VDW GmbH, Germany directly.

Problem	Possible Cause	Solution
THE DEVICE IS NOT FUNCTIONING PROPERLY	<ul style="list-style-type: none"> • Battery empty. • The charger is not correctly connected to the power socket. • The mains voltage being used to recharge the battery does not correspond to the voltage indicated on the adapter label. 	<ul style="list-style-type: none"> • Charge the battery. • Check whether the charger is plugged in correctly. • Check whether the original adapter is being used. • Reload factory default parameters (see details in Chapter 7.6 Default Parameters).
THE DISPLAY SHOWS INCORRECT VALUES	Unstable or fading display due to weak battery power.	Charge the battery.
THE MICROMOTOR WILL NOT START	Either the micromotor is improperly connected to the housing, or the VDW 6:1 contra-angle is defective.	<ul style="list-style-type: none"> • Check whether the micromotor cable is correctly connected to the micromotor housing. • Check whether the contraangle is working correctly. • Remove the contra-angle and set the maximum rotational speed; then restart the micromotor. • Calibrate without the contra-angle; then put the contra-angle back on the micromotor and repeat the calibration process.
THE FOOT PEDAL DOES NOT START THE MICROMOTOR	The foot pedal is either damaged or incorrectly connected.	First check whether the foot pedal is correctly connected to the device. Press the foot pedal again. If the micromotor does not start, then start the micromotor by pressing the ✓ Confirm key for 1.5 seconds. If this starts the micromotor, the foot pedal is defective. Contact your service centre to have the foot pedal replaced.

Problem	Possible Cause	Solution
THE BATTERY IS NOT WORKING PROPERLY	<ul style="list-style-type: none"> • Despite observance of all precautionary measures, the battery discharges too quickly. • The device only works if the charger is plugged into the power supply – however not in battery mode. 	The battery could be damaged. Send the device to your service centre.
CALIBRATION ERROR 1	An incorrectly-connected micromotor could disrupt the calibration process.	Check whether the micromotor is correctly connected.
CALIBRATION ERROR 2	A VDW 6:1 contra-angle with excessive resistance can cause the calibration process to abort.	<ul style="list-style-type: none"> • Check whether the contra-angle is malfunctioning. • Determine whether your contra-angle care and maintenance routines are correct. • To rule out the possibility of the motor being defective, recalibrate the motor without the contra-angle.
ROTARY FILE BLOCKED IN CANAL	<ul style="list-style-type: none"> • Incorrect file settings. • Too much pressure on the instrument. 	Switch to “backwards” ASR mode (red LED), restart the motor and carefully pull the file out.
RECIPROCATING FILE BLOCKED IN CANAL	<ul style="list-style-type: none"> • Too much pressure on the instrument. • File not cleaned often enough (see details in the RECIPROC® user’s manual). 	Remove the file using a pair of pliers by pulling it and carefully twisting it clockwise.
The display reads APEX LOCATOR ERROR 3	Apex locator cannot be switched on due to a malfunction.	<ul style="list-style-type: none"> • Check whether the lip clip and/or file clip cables are inserted incorrectly or damaged. • Check whether you have accidentally connected the file in the contra-angle to the file clip instead of the lip clip. • Check whether the original VDW 6:1 contra-angle is being used. • Check whether the original lip clip cable and the original charger with ferrite ring are being used.

Problem	Possible Cause	Solution
<p>APEX STATUS LED does not switch to:</p> <p>red (stays off):</p> <p>green:</p>	<p>Faulty connection between the lip clip cable and the unit. Defective lip clip cable.</p> <p>Lip clip and contra-angle file or lip clip and separate file clip connected too fast or incorrectly.</p>	<ul style="list-style-type: none"> • Contra-angle not functioning correctly. • Check whether the lip clip cable is inserted incorrectly or whether it is twisted or defective. • Clean the lip clip. • Check whether the measuring file has been correctly inserted into the contra-angle. • Check whether the contra-angle works. • Repeat the activation process and wait until the message “APEX LOCATOR ON” appears. • Alternatively, repeat the activation with the lip clip and file in the file clip, and wait until the message “APEX LOCATOR ON” appears.
<p>APEX LOCATOR IS OVERSENSITIVE – i.e., it displays the apex too quickly or indicates too short a working length.</p>	<ul style="list-style-type: none"> • Short-circuit due to excess fluid in the pulp chamber (rinsing solution, saliva, blood). • Direct contact between file and gingiva/ gingival proliferation, i.e. in cases with fractured metal crowns. • Direct contact between file and metal restoration work (crowns, parapulpal pins, amalgam fillings). • Lateral root canal. • Immature canal with large apex. 	<ul style="list-style-type: none"> • Dry the access cavity using an air blower or a cotton ball. • If heavy bleeding is present, wait until bleeding has been successfully stopped. • For isolation: <ul style="list-style-type: none"> - Adequate tooth build-up - Electrocauterisation • Insert a dental dam. • Use the protective silicon cover on the contra-angle. • Carefully widen the access cavity; add flow composites for isolation purposes if necessary. • Repeat length determination process. • Precise results not possible.

Problem	Possible Cause	Solution
<p>LENGTH DETERMINATION DOES NOT WORK, as the electrical measurement circuit cannot be closed.</p>	<ul style="list-style-type: none"> • Calcified or obliterated root canal. • Extremely dry root canal. • Old root filling residues creating a blockage/blockage through insufficient removal of intracanal medicament (i.e. calcium hydroxide). • As separate length determination: File too narrow for large root canal. • Lip clip incorrectly hung in patient's mouth. • As separate length determination: File incorrectly connected to the file clip – i.e. actually connected to the metal shaft. • Defective connection cable. 	<ul style="list-style-type: none"> • Check comparative X-rays for more information. If necessary, catheterize to working length with ISO 06/08 file. • Rinse with NaOCl solution; dry access cavity with cotton ball/air blower. • Do a comparative X-ray, and remove old gutta-percha/ intracanal medicament completely. • Length determination/ remove intracanal medicament completely. • If there is no contact to the file clip, select a larger file. • Place the lip clip in the patient's mouth again. • Check cable and socket connections again. • As separate length determination: check whether the file and the file clip are properly in contact with one another. • Clean the file clip using ethanol. • Check cables and sockets for visible damage.

Problem	Possible Cause	Solution
THE APEX LOCATOR WILL NOT TURN ON	<ul style="list-style-type: none"> • One of the connection cables could be inserted improperly or defective. • Connection error between the cables and the unit. • Contact between the lip clip and either the contra-angle file or the file in the file clip either too short or faulty. 	<ul style="list-style-type: none"> • Check all plug connections and retry. • If you tried to activate the apex locator by connecting the lip clip to the file in the contra-angle, try again using the lip clip and the file in the file clip. If this successfully activates the apex locator, then the connection between the lip clip and the file in the contra-angle is faulty. • Check the contact between the file and the lip clip. • If the results are still negative, this indicates that the contra-angle (only the original VDW 6:1 contra-angle may be used) is not transmitting signals correctly. The contra-angle needs to be replaced. • If the error persists following contact between the lip clip and the angle-piece file as well as between the lip clip and the file clip file, this could indicate a defective lip clip cable or a break in the measurement current circuit (file correctly inserted, file with electricity-conducting steel shaft). If the error cannot be fixed, switching the apex locator on and thus using it will not be possible.

11. Warranty

In addition to the warranty specified in the sales contract with the dental dealer, VDW GmbH provides the following factory service guarantee to customers directly:

1. VDW GmbH vouches for the correct construction of the product, the use of top-quality construction materials, the performance of all necessary testing and that the product is in compliance with all relevant laws and regulations.

The full functionality of the VDW.GOLD® RECIPROC® is guaranteed under a **36 month** warranty, beginning on the date the device is delivered to the customer (as indicated in the delivery papers containing the serial number of the product, provided by the supplier at the time of sale). The VDW 6:1 contra-angle has a **12-month** warranty.

The following wearing parts are not included in the warranty: 2-part measurement cable, lip clip and file clip.

The customer only has the right to guaranteed services during the warranty period, and only on condition that the VDW GmbH is informed of the defect in writing within a two month period from the date of discovery of the defect.

2. In case of a substantiated claim, the VDW Service Centre Munich will complete the repairs within 3 working days from the time the device arrives at VDW GmbH in Munich, not including time for return transport of the device to the customer.

3. This warranty covers only the exchange or repair of those individual components or parts which are affected by manufacturing defects. Costs associated with customer visits by dental-supply technicians will not be reimbursed by VDW GmbH, nor will packing costs incurred by the customer. The customer has no right to additional compensation by VDW GmbH which goes beyond repairs (such as damage compensation). This warranty does not cover compensation for direct or indirect personal injury or property damage of any kind. The customer has no right to compensation for damages due to device downtime.

4. The warranty does not extend to damage which VDW GmbH can prove is the result of user neglect in terms of normal device care (see user's manual), especially as regards battery charging, discharging and care in accordance with the user's manual instructions as well as careful and regular care of the VDW 6:1 contra-angle as discussed in the special VDW 6:1 contra-angle user's manual. The warranty expressly excludes defects which:

- occur as a result of damage during transport to VDW GmbH for repair purposes;
- caused by atmospheric events such as lightning, fire and humidity. This warranty is automatically invalidated if the product is inappropriately repaired, modified or otherwise handled by the user or by unauthorized third parties.

5. The warranty is only valid if the invoice confirming product delivery date is included with the device being sent for repair work.

6. Legal claims, for example those related to product liability rights, or claims against those from whom the customer obtained the product, especially against dental suppliers, remain unaffected.

12. Dr's Choice Tables

Record the file sizes and the appropriate values for your custom torque and rotation settings in the following table (For details see Chapter 7.5.2 Dr's Choice (only in rotary mode)):

File Position	File Type	gcm	rpm
01			
02			
03			
04			
05			
06			
07			
08			
09			
10			
11			
12			
13			
14			
15			

Preset Values in Dr's Choice Program on Delivery:

File Position	gcm	rpm
01	30	300
02	50	300
03	70	300
04	100	300
05	120	300
06	150	300
07	170	300
08	200	300
09	220	300
10	250	300
11	270	300
12	300	300
13	320	300
14	350	300
15	400	300

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Appendix

Electromagnetic Emissions and Immunity (English)	177
--	-----

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Electromagnetic Emissions and Immunity

The device is intended for use in the electromagnetic environment specified below. The user of the device should assure that it is used in such an environment.

Emissions		
Emission test	Conformity	Electromagnetic Environment - guidance
RF Emissions Cispr 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions Cispr 11	Class A	The device is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the device or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A Conforms	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Conforms	

Immunity Aspects			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	Test level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Burst/Fast Transient EN 61000-4-4	± 2 kV power supply lines	± 2 kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% U_T (> 95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 seconds	< 5% U_T (> 95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be that of a typical commercial or hospital environment.
Note: U_T is the mains voltage in AC before applying test level			

Immunity Aspects at r.f.

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the navigator should assure that it is used in such an electromagnetic environment.

Immunity test	Test level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
RF conducted EN 61000-4-6	3 Veff from 150 kHz to 80 MHz	3 Veff from 150 kHz to 80 MHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2 \cdot \sqrt{P}$ from 150 kHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$ from 80 MHz to 800 MHz $d = 2.3 \cdot \sqrt{P}$ from 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^{a)} should be less than the compliance level in each frequency range^{b)}. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
RF radiated EN 61000-4-3	3 Veff from 80 MHz to 2.5 GHz	3 Veff from 80 MHz to 2.5 GHz	

Notes:

- (1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ME Equipment is used exceeds the applicable RF compliance level above, ME Equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ME Equipment.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150 kHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$	From 80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$	From 800 MHz to 2 GHz $d = 2.3 \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

VDW.GOLD® RECIPROC®

CE
0123

Manufacturer
VDW GmbH • Bayerwaldstr. 15 • 81737 Munich • Germany
Phone +49 89 62734-0 • Fax +49 89 62734-304
info@vdw-dental.com • www.vdw-dental.com