Vacuum Pump VP5

EN Operating Instructions
Page 3

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ivoclar
vivadent[®]

Vacuum Pump VP5

	Page
List of Parts	4
1. Introduction / Signs and Symbols	5
1.1 Preface	
1.2 Introduction	
1.3 Signs and symbols	
2. Safety First	6
2.1 Indications	
2.2 Health and safety instructions	
3. Product Description	7
3.1 Components	
3.2 Functional description	
3.3 Hazardous areas and safety equipment	
4. Installation	8
4.1 Unpacking and checking the contents	
4.2 Selecting the location	
4.3 Making the connections	
5. Start-Up	9
6. Operation	9
7. Maintenance, Cleaning, and Diagnosis	10
7.1 Monitoring and maintenance	
7.2 Cleaning	
7.3 Notes on maintenance	
8. What If...	11
8.1 Technical malfunctions	
8.2 Repairs	
9. Product Specifications	12
9.1 Delivery from	
9.2 Technical data	
9.3 Acceptable operating conditions	
9.4 Acceptable transportation conditions	

List of Parts

- 1 Cable connection with bend protection
- 2 Power cord
- 3 Air vents
- 4 Vacuum hose connection
- 5 Housing
- 6 Vacuum hose
- 7 Transportation grip
- 8 Sound absorber
- 9 Sound absorber connection
- 10 Protective cap
- 11 Rating plate
- 12 Rubber feet

1. Introduction / Signs and Symbols

1.1 Preface

Dear Customer

Thank you for having purchased the Programat® Vacuum Pump. This apparatus is a high-quality, technical product. Please read the Operating Instructions carefully and use the apparatus according to the Instructions. Should you have any further questions, please do not hesitate to contact your local distributor or Ivoclar Vivadent AG directly.

1.2 Introduction

Apparatus: Vacuum Pump VP5
Target group: Dental lab personnel

This vacuum pump is suitable for evacuating the firing chamber of the Programat ceramic furnace. It is a quietly-running double-diaphragm pump with high vacuum performance.

These Operating Instructions are designed to assure the correct, safe, and economic use of the VP5 vacuum pump. They are divided into several, clearly structured chapters to help you find specific topics quickly and easily.

Signs and symbols (pictographs) are used to inform you quickly and clearly about health risks, important details, and contraindications.

We recommend keeping the Instructions in a safe place near the pump so that they are immediately accessible if necessary. Should you lose the Operating Instructions, extra copies can be ordered at a nominal fee from your local Ivoclar Vivadent Service Center or downloaded from www.ivoclarvivadent.com/downloadcenter.

1.3 Signs and symbols

The signs and symbols in these Operating Instructions and on the pump facilitate the finding of important points and have the following meanings:

Operating Instructions:

Risks and dangers

Information

Contraindication

Pump:

Alternating current

2. Safety First

This chapter is especially important for personnel who work with the VP5 vacuum pump or who have to carry out maintenance or repair work. This chapter must be read and the corresponding instructions followed.

2.1 Indications

This vacuum pump has been developed especially for the evacuation of Programat ceramic furnaces from Ivoclar Vivadent and it should be used for this purpose only.

Other uses than the ones stipulated are contraindicated. The manufacturer does not assume any liability for damage resulting from misuse. The user is solely responsible for any risk resulting from failure to observe these instructions.

Further instructions to assure proper use of the pump:


- The instructions, regulations, and notes in the vacuum pump's Operating Instructions must be observed.
- The pump must be operated under the indicated environmental and operating conditions (Chapter 9).
- The VP5 must be properly maintained.

2.1.1

Risks and dangers

Foreign objects or liquids must not enter the air vents, since this may block the motor or result in an electrical shock.

2.2 Health and Safety Instructions

- Make sure that the pump is only used for the indications stipulated in these Operating Instructions.
- Do not use the pump in areas where there is an explosion hazard.
- Components connected with the pump must comply with the pneumatic data of the pump.
- Observe the corresponding safety regulations when connecting the pump to the power supply.
- If the operation of the pump is interrupted by the thermostat because of over-heating, the pump starts automatically once it has cooled down. Make sure that no danger results from such an incidence.
- Use only original spare parts.
- The pump complies with the following directives:
 - Directive 2006/42/EC on machinery
 - Directive 2004/108/EC on electromagnetic compatibility
- If the pump has been stored at very low temperatures or high atmospheric humidity, the pump has to be dried or left to adjust to the room temperature for approx. 1 hour (do not connect to the power yet).

Do not work with liquids near the pump. Should a liquid accidentally enter the pump, disconnect the power and consult Customer Service. Do not operate the pump.

- The pump is tested for use at altitudes of up to 2000 m above sea level.
- The pump may only be used indoors.
- The unit must not be technically modified.

3. Product Description

Any disruption of the protective conductor either inside or outside of the pump or any loosening of the protective conductor connection may lead to danger for the user in case of malfunction. Deliberate interruptions are not tolerated.

Do not use furnaces of other manufacturers that are not coordinated with the pump. Ask the corresponding manufacturer for further details. The vacuum pump has been especially developed for use with the Programat ceramic furnaces and must be used for this purpose only.

Disposal:

The pump must not be disposed of in the normal domestic waste. Please correctly dispose of old furnaces according to the corresponding EU council directive. Information regarding disposal may also be found on the respective national Ivoclar Vivadent homepage.

3.1 Components

The VP5 vacuum pump consists of the following components:

- Vacuum pump with power cord
- Vacuum hose (accessory)

3.2 Functional description

The diaphragm squeezed between the cylinder head and the housing is moved up and down by a con-rod. By means of this movement and by suction and eduction valves in the cylinder head, the firing and/or press chamber of the furnace is evacuated and the vacuum is built up.

3.3 Hazardous areas and safety equipment

Description of the hazardous areas of the pump:

Hazardous area	Type of risk
Fuse holder	Electrical shock
Ventilator	Injury

Safety equipment	Protective effect
Protective conductor	Protection from electrical shock
IP20 (protection type)	Protection from electrical shock
Thermoswitch	Protection from overheating
Ventilator cover	Protection from electrical shock and injury

4. Installation

4.1 Unpacking and checking the contents

Remove the pump from its packaging and check it for any transportation damage. We recommend keeping the original packaging for any transportation purposes. Use only the original packaging for shipping.

The electrical installation in the room where the pump is used must comply with the local and IEC requirements.

4.2 Selecting the location

Place the pump on a flat surface using the rubber feet (12). Make sure the pump is not exposed to direct sunlight. Do not place the pump in the immediate vicinity of heaters or other sources of heat. Make sure that there is enough space for air circulation between the wall and the air vents (3).

Make sure to select the location in such a way that the sound absorber is not knocked off by accidentally abutting against laboratory furniture.


4.3 Making the connections

Checking the rating plate and power supply

Make sure that the voltage indicated on the rating plate (11) complies with the local power supply. Should this not be the case, do not connect the pump.

Connecting the vacuum hose

Remove the protective cap (10) from the vacuum hose connection (4) and slip the vacuum hose on the vacuum hose connection (4).

Connect the other end of the vacuum hose with the furnace.

Please check if the hose is properly connected.

Connecting the power

Plug the power plug of the power cord (2) into the corresponding socket in the Programat. This is the power supply at the same time.

5. Start-Up

The vacuum pump is automatically controlled by the furnace, once the power cord of the pump has been connected with the Programat furnace.

6. Operation

The pump is operated by the furnace programs.

7. Maintenance and Cleaning

For safety reasons, disconnect power plug before maintenance and cleaning.

7.1 Monitoring and maintenance

Basically, the diaphragm vacuum pump does not require maintenance unless the performance deteriorates due to contamination. Therefore, conduct a vacuum test at least once a year.

The time for maintenance procedures depends on the frequency of use and the location of the pump. For that reason the recommended times are only approximates.

What	Part	When
Check all plug-in connections for correct fit.	Connections	Weekly
Check the hose to make sure that no condensation water has accumulated in the system.	Vacuum hose	Weekly

This apparatus has been developed for typical use in the dental laboratory.

If the product is used in a production facility, for industrial applications, or in continuous operation, premature ageing of certain spare parts have to be expected.

These spare parts are e.g.:

- Sound absorber
- Diaphragm
- Valve reed

These spare parts are not covered by the warranty.

Please also observe the shorter service and maintenance intervals.

7.2 Cleaning

Remove dust from the pump from time to time. Use a vacuum cleaner with a cleaning brush for that purpose.

Item	Frequency	Cleaning material
Housing	if required	soft, dry cloth
Sound absorber (8)	if dirty	sound absorber has to be replaced

7.3 Notes on maintenance

The diaphragm may only be changed by an Ivoclar Vivadent Service Center. The fuse must be changed by qualified personnel.

Please observe the relevant rules and regulations for the prevention of accidents as well as the safety and medical regulations.

8. What If ...?

8.1 Technical malfunctions

Malfunction	Double-check	Action
1 The pump is not working	<ul style="list-style-type: none">– Is the pump correctly connected?– Has the thermostat been activated because of overheating?– Sound absorber clogged	<ul style="list-style-type: none">– Connect pump correctly.– Allow the pump to cool and switch on again.– Check pump location for adequate circulation of air.– If the problem cannot be solved, contact your local Ivoclar Vivadent Service Center.– Replace sound absorber.
2 The vacuum performance is declining or is bad	<ul style="list-style-type: none">– Has the test program for the pump been carried out?	<ul style="list-style-type: none">– Carry out the test program for the pump (PX1).– Send vacuum pump to Service Center for maintenance.
3 Vacuum pump is unusually noisy	<ul style="list-style-type: none">– Is the sound absorber defective or incorrectly mounted?– Is there a leak in the furnace system or the pump?	<ul style="list-style-type: none">– Connect sound absorber correctly or replace it, if it is defective.– Check connections and hose for tightness.– Clean sealing of the furnace head.– Check pump performance with the corresponding test program.

8.2 Repairs

Repairs may only be carried out by the personnel of a certified Ivoclar Vivadent Service Center. Please contact your Customer Service representative.

9. Product Specifications

9.1 Delivery form

- Vacuum pump VP5
- Operating Instructions
- Warranty Card
- Declaration of conformity

9.2 Technical data

Power supply:	230 V / 50–60 Hz / 0.40 A 200 V / 50–60 Hz / 0.50 A 115 V / 50–60 Hz / 0.85 A 100 V / 50–60 Hz / 0.95 A
Tolerated voltage fluctuations:	+/- 10%
Power consumption:	80 W
Suction strength:	15 l / min.
Final vacuum (measured with a 4-liter test container):	approx. 20 mbar
Dimensions:	415 x 122 x 230 mm (length, width, height)
Weight:	8.0 kg
Types of fuses:	230 V: T 1.0 A 250 V 200 V: T 1.25 A 250 V 115 V: T 2.5 A 250 V 100 V: T 2,5 A 250 V

Safety information :

- The pump complies with the following directives:
- Directive 2006/42/EC on machinery
 - Directive 2004/108/EC on electromagnetic compatibility

9.3 Acceptable operating conditions

- Ambient temperature range of:
+5 °C to +40 °C (+41 °F to +104 °F)
- Acceptable height:
The pump is tested for use at altitudes of up to
2000 m above sea level.
- Pollution degree 2 Overvoltage category II

9.4 Acceptable transportation conditions

- Ambient temperature range of:
-10 °C to +55 °C (-4 °F to +131 °F)
- Maximum relative humidity: 80%

This apparatus has been developed for use in dentistry. Setup and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is solely responsible for testing the apparatus for its suitability for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.