

Contents

1. Introduction	04
- Features - Application	
2. Safety information	05
SymbolsNotes and referenceInstallation and connectionOperationMaintenance	
3. List of components	09
4. Parts and functions	10
5. Installation and connection	11
6. Direction for use	13
7. Selection of curing modes	14
8. Maintenance	15
Cleaning and sterilizationMaintenance and replacement of the Batte	ery
9. Trouble shooting	16
10. Specification	17
11. Accessories	18
12. Warranty	18
13. Warranty card	19

1 Introduction

Thank you for purchasing the device.


The instructions described below have been designed to provide necessary information for successful operation, daily maintenance and safety for use of this unit.

Features

- Dual wavelength LED (peak 405nm & 460nm)
- Autoclaveable guide-cap (134 °C)
- Variety of curing modes
- Easy-to read OLED screen
- Enhanced capacity battery
- Lithium-ion battery
- Interchangeable lithium-ion battery
- Chargeable battery alone
- Automatic sleep mode
- Ergonomic and familiar design

Application

 Polymerization of light-cure catalyst containing dental materials with a wavelength spectrum of 400 ~ 490nm.



2 Safety information

It is important to understand the following information to use this product safely.

Symbols



Indicates a potentially hazardous situation which, if not avoided, may result in serious injury or death.



Indicates a potentially hazardous situation which, if not avoided, may result in minor/moderate injury or in damage of property

Notes and References

Notes and References used in this manual are shown below.



Indicates that potential malfunctions may occur if the description is ignored.



Supplementary information to facilitate better operation.

Installation and Connection



- The adaptor should be disconnected from the charger before connecting the power cable to prevent an electric shock.
- Keep away from sources of ignition to avoid fire or explosion.
- Protect the handpiece, charger, and battery from any liquid to avoid electric shock or functional error.
- Place the charger on a flat and stable surface. Do not drop handpiece or charger, in order to avoid internal defect that may cause overheating or fire.
- Do not place the electrical cord near a gas burner other sources of heat.
 Should the cord catch on fire, replace with a new one to avoid short-circuit.



- Do not use or leave the unit in areas of high temperature or exposed to direct sunlight;
 internal circuit defect may cause overheating or fire.
- Securely connect the dust-free power plug to the outlet.

 Use this device in well ventilated areas. Overheating may cause fire.
- Do not touch the terminal areas of charger or battery with fingers or conductive objects, such as metal wire or safety-pin. It may cause electric shock or fire.
- Do not use any charger other than supplied. It may result in electric shock or fire.

2 Safety information

Operation



- Do not use the device without guide-cap. It may result in burn cause by heat from LED.
- Do not plug the power cord with wet hands. This may result in an electric shock.
- In case of fluid leak from battery or partial discoloration/deformation of the Handpiece, immediately discontinue use.
- Place the battery in correct direction when charging. Wrong directed placement may cause malfunction, explosion, or fire.
- In case of partial discoloration of the charger, immediately disconnect the power plug and repair the charger.
- In case battery fluid contacts the eye, immediately flush with copious amount of clean water and seek medical attention.
- Do not use this product on patients with a history of cataract, retinal damage or allergic reaction to light.
- Do not use this product on a person with a pacemaker as it may cause malfunction of the pacemaker.
- Avoid exposing the light on intraoral soft tissue as it may affect the soft tissue.
- Use carefully the guide-cap to avoid broken lens on the guide-cap.
- Do not put the battery on the charger upside down. This may cause an explosion or fire.



- Disconnect the power cord by grasping the power plug to avoid injuries or burns.
- Avoid damage or use of unnecessary force on the power cord. This may cause an electric shock or fire.
- Should battery fluid contact the skin or clothes, immediately wash it off with copious amounts of water to avoid skin damage.
- Take the battery out of the handpice if it won't be used for more than 1 week.
- Connect the power plug only when battery is being charged. Insulation degradation may cause electric shock or fire.
- Do not use this product for any purposes other than dental use.
- The product may malfunction or turn off due to static electricity.

2 Safety information

Maintenance



- Use only the battery supplied to avoid fluid leak or explosion.
- Never disassemble or overhaul the product. Unusual performance may cause injuries, electric shock, or fire.
- Ensure that charger and battery are always separated from the handpiece prior to cleaning.



- If the OLED screen malfunctions, remove the battery from the handpiece and reinsert it.
- User must be in charge of operation, maintenance and check.
- Take care a lens on the guide-cap to prevent broken when cleaning.


3 List of components

User manual

Please check the contents described below when removing the contents from the original packaging.



4 Parts and functions



① Handpiece

Start button : Operation and turn On/Off
 Mode button : Curing mode select
 Control button : Curing time select
 OLED screen : Display device state

② Guide-cap: Autoclaveable, easy to replace

3 Battery

4 Charger: Recharge battery display charging process

⑤ Adapter


6 Power cable

5 Installation and connection

Refer to the information supplied on safety precautions during installation and connection.

1. Connection between adapter and charger

- Securely insert the power cable with 100~240V power outlet.
- 2) Connect power cable with adapter.
- 3) Connect adapter with charger





- Use only the original adapter supplied.
- When the adapter is connected with charger, red and blue light will blink.
- If red and blue light is not blinking when adapter is connected, it is a charger malfunction.

2. Inserting battery


- ① Carefully insert the battery into the handpiece as directed by the arrow.
- ② The battery turns with the arrow, as seen in image. *Avoid forcing the battery.



5 Installation and connection

3. Inserting guide-cap

 Slowly insert the guide-cap into the handpiece as picture. Insert it till right position.





- Replace the guide-cap when it is dirty, or intensity is lower than normal condition cause by scratch or broken on lens.
- For convenience, sterilize the guide-cap before use.

4. Charging

- Insert the handpiece into the vertical slot at the front of the charger.
- Red light on : Charging in progress
- Blue light on : Charging process is complete





 If the remaining battery is extremely low, red light may flicker while recharging, which is not a malfunction.




- It will take approximately 3~4 hours before the battery will be fully charged.

6 Direction for use

1. Activating handpiece

 Press any button, then OLED screen will be on standby mode and time display used last.



Select curing mode and time setting with mode button and control button.



Initial screen



Standby screen (last used mode & time)





- When battery inserted, the unit is automatically turned on.

2. Start of exposure

Press start button to start exposure with the selected mode.



- Press any button if you wish to end the curing process.
- Alarm during operating informs time passage.

3. Battery icon

This icon indicates the power remaining in the battery. Recharge or replace the battery when "LOW BATTERY" sign appears on the display screen.



 Sleep mode: If unit is left unused for a few minutes, all functions are automatically transferred to sleep mode with a short beep.
 Hold down start button for approximately 3-5 seconds to shut off light. If you want to turn on, press any button.

- Set-up memory: The handpiece will record the setting used last.


7 Selection of curing mode

Curing mode can be set by pressing the mode button. Various curing programs are available in this unit.

1. High mode


Time setting: 5, 10 sec.

Quick polymerizing with high intensity



2. Low mode


Time setting: 10, 20, 30 sec. Polymerizing with low intensity



3. Soft start mode

Time setting: 20 sec. (fixed)


Special mode to avoid polymerization shrinkage



4. Pulse mode

Time setting: 10, 20 sec.

Special mode to avoid polymerization shrinkage



Cleaning and sterilization

- Clean the handpiece and charger gently with clean, soft, lint-free cloth.
- Avoid using solvents which contain ethanol or acetone.
- Do not use sharp or pointed tools to clean as they may damage or scratch the plastic.
- Clean the guide-cap with a soft cloth and disinfect. The guide-cap is autoclaveable.
- Remove restorative materials stuck to the guide-cap gently with a plastic spatula.
 Avoid using sharp tool as they may scratch the surface of the lens and affect the intensity of light.

Maintenance and replacement of the battery

- Exposure to high room temperature may shorten the battery life and affect device performance.
 Use the device in normal temperature.
- It is recommended for rechargeable battery to be charged immediately before use to maximize efficiency.
- Remove the battery from the handpiece, if not being used for a long time.
- If the running time declines with a full-charged battery, replace the battery before use.

Trouble shooting

If your curing light unit is not operating normally, please refer to the following points.

	Trouble	Cause and solution	
1	Control button does not work	If the control button is pressed with excess force, it will not work. In such case press the button gently.	
2	Display screen does not work	If the battery is completely discharged or connected incorrectly, display will not appear. Recharge the Battery or reconnect the Battery to reset the device.	
3	Power off while in operation	If the start button is pressed for a moment, power will be turned off. Press any button again. Check the battery connection again. If the battery is completely discharged, device can be turned off. Recharge the battery.	
4	Not able to charge	If the power cord or etc. is not correctly connected, the battery can not be charged. If the battery or handpiece is not properly placed on the charger, the battery can not be charged.	
5	Poor intensity	Dirty, damage of guide-cap can make poor intensity. Clean the guide-cap or replace new one.	

^{*} If the trouble is not settled after upper solutions, please contact local dealer or Good Doctors Co., Ltd.

10 Specification

	Operation voltage	3.7V						
Handpiece	Peak wavelength	405nm & 460nm						
piece	Dimensions (included battery)	Approx. 232 x 40 x 26mm						
	Weight (included battery)	Approx. 107g						
0	Operation voltage	5V						
Charger	Dimensions	Approx. 110 x 105 x 68mm						
ïr	Weight	Approx. 100g						
	Туре	Lithium Ion battery						
Battery	Nominal voltage	DC 3.7V						
tery	Dimensions	Approx. 75 x 25 x 25mm						
	Weight	Approx. 52g						
Guide-cap	Dimensions	Approx. 64 x 16 x 16mm						
э-сар	Weight	Approx. 5g						

■ Operation, transport, and storage condition for the main unit

Input voltage

Output voltage

	Operation	Storage	Transport
Temperature	10~40C	0~40C	-10C~50C
Humidity	10~80%	0~80%	0~80%
Atmospheric pressure	800 ~ 1060 hPa		

^{*} This specification is subject to change without prior notice

AC 100 ~ 240V

DC 5V

Accessories





Guide-cap (x 3)

Battery

12 Warranty

This product is shipped after strict inspection. Should the product have any problem under normal usage within one year from the date of purchase, repair will be free of charge.

However, repair of consumables (guide-cap, battery, etc.) are not covered by the warranty. For more information, refer to the warranty card attached.

Serial number check.





EC REP Good Doctors Germany GmbH

Gerhard-Domagk-Str.2, 53121 Bonn, Germany TEL: +49-228/53441465 FAX: +49-228-53441466 Website: www.gooddrs.de

Good Doctors Co., Ltd.

#208, B-dong, 283 Bupyeong-daero, Bupyeong-gu (Woolim Lions Valley ,Cheongchoen-dong), Incheon, 21315. Republic of Korea

Tel: + 82 32 424 6325 Fax: + 82 32 424 6326

E-mail: good@gooddrs.com Website: www.gooddrs.com

. उटवर्नेठाः ५ उटवर्नेठाः ५ उटवर्नेठाः ६ उटवर्नेठाः ६ उटवर्नेठाः ६ उटवर्नेठाः 50000



We guarantee the quality of the product as follows

Name of product:

Name of model:

1. If there is any defect of manufacturer or product has any fault itself within the assurance period. it will be replaced with the same product.

< Assurance details >

Customer information

- Office phone: - Office name :

- Dealer name : Distributor

 Malfunction and damage according to user's intention or accident. Malfunction and damage according to using incompaible product The manufacturer does not accept any liability resulting misuse.

- Malfunction and damage according to natural disaster such - Malfunction and damage according to user's carelessness. 2. Below cases are not effective for the free refund.

as fire, earthquake and flood etc.

Assurance period

- Date of purchase : from

Product is for 1 year except accessories.

, 50

Good Doctors Co., Ltd.

#208, B-Dong, 283 Bupyeong-daero, Bupyeong-gu, Incheon, Korea (Woolim Lions Valley, Cheongcheon-dong) (Zip code : 21315) Tel: +82 32 424 6325 Fax: +82 32 424 6326

وموظفرت وموظفرت وموظفرت وموظفرت وموظفرت وموظفرت وموظفرت

