mectron

medical technology

EN

IT

DE

FR

ES

SV

- Use and maintenance manual
- Manuale d'uso e manutenzione
- Gebrauchs- und Wartungshandbuch
- Mode d'emploi et maintenance
- Manual de uso y mantenimiento
- Bruks-och underhållsanvisning

PIEZOSURGERY[®] touch

CE
0476

→ *Fold out this page*

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



→ *Ouvrir la page*

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→ *Vik ut denna sida*

→ Use and maintenance manual

PIEZOSURGERY[®] *touch*

→ TABLE OF CONTENTS

01 →	INTRODUCTION	4
01.1	Intended use of the PIEZOSURGERY® <i>touch</i>	4
01.2	Description of the device	5
01.3	Responsibility waiver	5
01.4	Safety precautions	6
01.5	Symbols	8
02 →	IDENTIFICATION DATA	9
02.1	Identification plate of the device	9
02.2	Identification data of the handpiece	9
02.3	Identification data of the inserts	9
03 →	DELIVERY	10
03.1	List of the components of the PIEZOSURGERY® <i>touch</i>	10
04 →	INSTALLATION	11
04.1	First installation	11
04.2	Safety precautions during the installation	11
04.3	Connection of the accessories	12
05 →	USE	14
05.1	Switching the device on and off	14
05.2	Description of the keyboard	15
05.3	Safety precautions before and during use	18
05.4	Instructions for use	19
05.5	Important information on the inserts	22

06	→ MAINTENANCE	23
07	→ DISPOSAL MODES AND PRECAUTIONS.....	23
08	→ TECHNICAL DATA	24
	08.1 Electromagnetic compatibility	25
09	→ TROUBLESHOOTING	29
	09.1 Diagnostic system and symbols on the keyboard.....	29
	09.2 Quick solution to problems	30
	09.3 Replacement of the fuses	32
	09.4 Sending the device to an Authorized Mectron Service Service ...	33
10	→ WARRANTY	34
11	→ BIBLIOGRAPHY	205

01 → INTRODUCTION

EN

Carefully read this manual before proceeding with the installation, use, maintenance, or other operations on the device. Always keep this manual within reach.

IMPORTANT: to prevent harm to persons or damages to objects, read all the "Safety precautions" present in the manual with special care. Depending on their degree of seriousness, the safety precautions are classified with the following indications:

-  **WARNING**
(always refers to personal injury)
-  **CAUTION**
(refers to possible damage to property)

The purpose of this manual is to make the operator knowledgeable of the safety

precautions, the installation procedures, and the instructions for a correct use and maintenance of the device and its accessories.

Use of this manual for purposes other than those strictly tied to the installation, use and maintenance of the device is forbidden. The information and illustrations in this manual are updated as of the date of issue reported on the last page.

Mectron is committed to continuously update its products with possible modifications to device components. In case you uncover discrepancies between what is described in this manual and the device in your possession, contact your Retailer or the After-Sales Service of Mectron for clarifications and support.

01.1 → INTENDED USE OF THE PIEZOSURGERY® *touch*

The PIEZOSURGERY® *touch* is a piezoelectric ultrasonic device intended for bone surgery in the following applications:

- osteotomy and osteoplasty techniques;
- implantology;
- periodontal surgery;
- orthodontic surgery;
- endodontic surgery;
- maxillofacial surgery

The PIEZOSURGERY® *touch*, and its inserts, can operate as a piezoelectric scaler in the following dental applications:

- scaling: all the procedures for the removal of bacterial plaque deposits and supra-gingival, sub-gingival, and interdental

calculi, and for the removal of stains;

- periodontology: periodontal treatment for the scaling and root-planing/debridement without damage to the periodontium, including cleaning and irrigation of the periodontal pocket;
- endodontics: all treatments for root canals preparation, irrigation, filling, condensation of the gutta-percha and retrograde preparation;
- restoration and prosthesis: cavity preparation, removal of denture, condensation of the amalgam, finishing of the cervical margin, inlay/onlay preparation.

 **WARNING:** The device must be used in a dental practice or ambulatory. The device cannot operate in environments where anesthetic or flammable mixtures are present.

 **WARNING: Qualified and specialized personnel.** The device must be used exclusively by specialized and appropriately trained personnel such as a Surgeon. The

use of the device does not cause side effects if it is used correctly. An improper use might cause tissues heating.

 **WARNING:** Use the device only for the intended use.

Failure to observe this precaution may cause serious injuries to the patient, the operator, and damages/breakdowns to the device.

01.2 → DESCRIPTION OF THE DEVICE

The user interface has been optimized with the PIEZOSURGERY® *touch* making all the functions readily available by integrating them in the touch keyboard.

The PIEZOSURGERY® *touch* is a device that uses ultrasonic piezoelectric technology to generate mechanical microvibrations of the inserts (from 20 to 60 µm), to effectively cut mineralized tissues. This allows an efficient and safe cutting which preserves the integrity of the osteotomized surfaces.

The micrometric, ultrasonic vibrations of the inserts provide greater precision and a selective cutting action compared to traditional methods such as drills or oscillating saws (which act with macrovibrations), therefore minimizing traumatic effect on soft tissues .

The cavitation effect of the irrigating solution helps to keep the operatory field blood-free. This provides an optimal intra-operative visual control thus increasing safety, even in areas that are anatomically most difficult to access.

01.3 → RESPONSIBILITY WAIVER

The manufacturer Mectron disclaim any liability, expressed or implied, and shall have no responsibility for any direct, indirect or other damages and personal injury arising out in connection with any practice in the use of the device and its accessories.

The manufacturer Mectron shall be under no liability, expressed or implied, with respect to any damages (personal injury and/or damage to property) which might arise or be caused, whether by the customer or by any of the users of the product and its accessories, as result of:

- 1 Use or procedures different than those specified in the intended use of the product;
- 2 The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 08 - TECHNICAL DATA;
- 3 The device is not used in compliance with all the instructions and precautions described in this manual;
- 4 The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions;
- 5 The assembly operations, extensions, adjustments, updates, and repairs on the device are performed by personnel not authorized by Mectron;
- 6 Improper use, mistreatments, and/or incorrect interventions;
- 7 Any and all attempts to tamper with or modify the device, under any circumstance;
- 8 Use of non-original Mectron inserts that entail a finite damage to the threading of the handpiece, thus compromising correct operation and causing risk of harm to the patient;
- 9 Use of non-original Mectron inserts, used in accordance to designed and tested settings of Mectron original inserts. The correct use of the settings is guaranteed only with original Mectron inserts;
- 10 Lack of stock materials (handpiece, inserts, wrenches) to be used in the event of device stop due to fault or of inconveniences.

01.4 → SAFETY PRECAUTIONS

EN

⚠ WARNING: Do not install the device in places where there is a risk of explosion.
The device cannot operate in environments where anesthetic or flammable mixtures are present.

ⓘ CAUTION: In the case that the end user, when operating in his or her own medical study or clinic, must subject the electro-medical equipment and systems to periodical inspections in order to adhere to imposed requirements, the test procedures that must be applied to electro-medical equipment and systems to evaluate safety must be performed in line with norm EN 62353 'Electro-medical devices - Periodical inspections and tests to be performed after repair interventions on electro-medical devices'.

⚠ WARNING: Checking device status before the treatment.
Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

ⓘ CAUTION: The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

ⓘ CAUTION: To avoid any risk of electric shock this device must be grounded.

⚠ WARNING: When the lid of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® *touch* must not be operated. The parts in movement may cause harm to the operator.

⚠ WARNING: All the accessories of the new or repaired devices are not sterile. All new and repaired accessories are supplied in non-sterile conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with

the instructions given in the Cleaning and Sterilization Manual.

⚠ WARNING: Infections control.
To ensure maximum safety of the patient and the operator, before using all the reusable parts and accessories, make sure you have first cleaned and sterilized them, following the instructions of the Cleaning and Sterilization Manual.

ⓘ CAUTION: Contraindications. Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ WARNING: Breakage and wear-out of the inserts. High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert.

Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used. If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them.

The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent the patient from ingesting fragments of broken inserts.

When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out. Use of a worn-out insert reduces the cutting performance and can cause necrosis of the bone surface treated. During the intervention, frequently check that the insert is intact, especially in its apical part.

During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use. Do not exert excessive pressure on the inserts during their use.

⚠ WARNING: Only use original Mectron inserts, accessories, and spare parts.

⚠ WARNING: Contraindications. Do not use the PIEZOSURGERY® *touch* on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ CAUTION: Contraindications. Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ WARNING: Contraindications.
Interference from other equipment. An electrical scalpel or other electro-surgical

units near the PIEZOSURGERY® *touch* device may interfere with its correct operation.

⚠ WARNING: Contraindications.
Interference with other equipment. Though compliant to the standard IEC 60601-1-2, the PIEZOSURGERY® *touch* may nonetheless interfere with other devices nearby. The PIEZOSURGERY® *touch* must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

⚠ CAUTION: No modification of this equipment is allowed.

01.5 → SYMBOLS

EN

Serial number

Lot number

Product code

CAUTION: read the instructions for use

Operating instructions

Temperature limitation - transport and storage conditions

Humidity limitation - transport and storage conditions

Atmospheric pressure limitation - transport and storage conditions

It indicates compliance with the CE 93/42 CEE EN 60601-1 and EN 60601-1-2.
Notified body: CERMET.

MET Mark
UL-CSA conformity

Manufacturer

Do not allow fingers to contact moving parts

QTY.1

Quantity of parts in the pack = 1

Single-use

Applied part of type "B" as per norm EN 60601-1

Can be sterilized in autoclave up to a maximum temperature of 135° C

Non-sterile

The device and its accessories must not be disposed of or treated as solid urban wastes

Biohazard

I

Activation switch "on"

0

Activation switch "off"

Alternating current

Connection of the foot pedal

Equipotentiality

For US market only

Rx Only

CAUTION US Federal Law restricts this device to sale by or on the order of a dentist

02 → IDENTIFICATION DATA


An exact description of the model and of the serial number of the device will enable our After-Sales Service to provide fast and efficient support.

Always refer these data whenever you contact an Authorized Mectron Service Center.

EN


02.1 → IDENTIFICATION PLATE OF THE DEVICE

Every device has an identification plate that bears the technical characteristics and the serial number. The identification plate is located on the inferior panel of the device. Additional specifications are reported in this manual (See Chapter 08 - TECHNICAL DATA).


02.2 → IDENTIFICATION DATA OF THE HANDPIECE

Each handpiece is laser-marked with its serial number (ref. 1), and Mectron and PIEZOSURGERY® *touch* logos (ref. 2).

02.3 → IDENTIFICATION DATA OF THE INSERTS

The following data are laser-marked on each insert: the name of the insert (ref. 3), the Mectron logo (ref. 4) and the lot number to which the insert belongs (ref. 5).

03 → DELIVERY

03.1 → LIST OF THE COMPONENTS OF THE PIEZOSURGERY® *touch* See inside cover

EN

PIEZOSURGERY® *touch* consists of:

- A device body
- B peristaltic pump
- C irrigation bag support rod
- D PIEZOSURGERY® *touch* handpiece fixed support
- E PIEZOSURGERY® *touch* handpiece mobile support
- F manual of use and maintenance and manual of cleaning and sterilization
- G electrical power supply cable
- H foot pedal with bracket, cord and connector
- I touch keyboard protective folio
- J handpiece complete with cord, LED light front cone and metal cone
- ⚠ **CAUTION** the handpiece and the cord cannot be separated
- K flow adjuster
- L inserts/insert kit
- M torque wrench
- N peristaltic pump tube
- O bone grafting kit (optional)
- P case

The PIEZOSURGERY® *touch* consists of accessories that can be ordered separately. Refer to the "Packing List" included in your package, so that you may exactly know the quantity and type of the accessories supplied with the device you have purchased. The package of the device is sensitive to strong collisions, because it contains electronic components. Therefore, special precautions must be taken for transport and storage. Do not overlap multiple boxes, in order not to squash the packages underneath. All the materials shipped by Mectron have been inspected upon their delivery. The device is delivered duly protected and packed.

When receiving the device, check for the possible presence of damages incurred during the transport and, should that be the case, file a complaint with the transporter. Preserve the package in the event that you need to send the device to an Authorized Mectron Service Center and so that you can place it back inside during prolonged periods of disuse.

⚠ WARNING: Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

04 → INSTALLATION

EN

04.1 → FIRST INSTALLATION

The device must be installed in a comfortable place suitable for its use.

The PIEZOSURGERY® *touch* can be purchased ready for it to be used instantly or may have to be enabled by typing an activation key.

In the case that your device requires an activation key, the procedures to follow may vary from country to country.

Always refer to your retailer to obtain information on the above.

04.2 → SAFETY PRECAUTIONS DURING INSTALLATION

⚠ WARNING: Contraindications. Interference with other equipment. Though compliant with the standard IEC 60601-1-2, the PIEZOSURGERY® *touch* may nonetheless interfere with other devices nearby.

The PIEZOSURGERY® *touch* must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

⚠ WARNING: Contraindications. Interference from other equipment. An electrical scalpel or other electro-surgical units near the PIEZOSURGERY® *touch* device may interfere with its correct operation.

⚠ CAUTION: The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

⚠ CAUTION: To avoid any risk of electric shock this device must be grounded.

⚠ WARNING: Do not install the device in places where there is a risk of explosion. The device cannot operate in environments where anesthetic or flammable mixtures are present.

⚠ CAUTION: Always position the device in way so that the power switch is easily reachable, since this switch is considered as a load-break switch.

⚠ WARNING: When the lid of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® *touch* must not be operated. The parts in movement may cause harm to the operator.

⚠ WARNING: Install the device in a place protected against collisions or against accidental sprays of water or liquids.

⚠ WARNING: Do not install the device above or near heat sources. Foresee adequate air circulation around the device when installing it. Leave adequate space, especially near the fan placed on the back part of the device.

⚠ CAUTION: Do not expose the device to direct sunlight or to sources of UV light.


⚠ CAUTION: The device can be transported, but it must be handled with care when it is displaced. Position the foot pedal on ground, so that it can only be activated intentionally by the operator.

⚠ CAUTION: Before connecting the handpiece cord to the device, make sure that the electrical contacts are perfectly dry. If need be, dry them with compressed air.

04.3 CONNECTION OF THE ACCESSORIES


Insert the irrigation bag support rod in the dedicated hole;

1


Insert the handpiece fixed support in the housing;
PLEASE NOTE The handpiece support can be placed in 4 different locations: lateral right, front right, front left, lateral left;

2


Take a peristaltic pump tube in silicon and insert it on the connector of the handpiece cord;

3


Insert the connector of the handpiece cord with the writing UP facing upward inside the connection in the front part of the PIEZOSURGERY® touch;


4

Insert the silicon tube in the peristaltic pump, proceeding as follows:
- open the pump cover completely;

5


6

- position the tube in the pump;
- close the pump cover completely;


⚠ WARNING: When the lid of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® *touch* must not be operated. The parts in movement may cause harm to the operator.

EN


7

- Connect the silicon tube end to the flow adjuster;

8

- Connect the flow adjuster to the liquid bag used for the treatment (bag not supplied). Hang the bag to its specific support rod;


9

- Connect the foot pedal to the back of the device in the socket marked with the symbol  by means of the plug of the pedal cable, until you hear a "click" sound;


PLEASE NOTE The foot pedal comes equipped with a bracket that allows it to be moved to the place most suitable for the operation, without the need to use your hands to move it;

The bracket can also be positioned horizontally if it is not used;

Insert the power supply cable in its connection located at the back of the device. Connect it to the power outlet;


10

Equipotential plug: The device is equipped with an additional equipotential plug located on the rear of the console. This plug is in accordance with DIN 42801. Insert the connector of the equipotential cord (optional) to the plug on the rear of the device's console.


The purpose of additional potential equalization is to reduce differences of potential which can occur during operation between the device's body and conductive parts of other objects within the medical environment;

11

Clean the touch keyboard with a soft cloth with low fiber-release and position the sterile protective foil on the surface of the clean touch keyboard.

12

05 → USE

05.1 → SWITCHING THE DEVICE ON AND OFF

Turning the device on


With the device in view from its front side, position the switch on the left of the device body on position "I", being careful not to press the foot pedal.

4 symbols appear on the device (ref. U inside cover) that turn off one by one. At this point, the device sets on the default setting and is ready for it to be used.


Turning the device off

With the device in view from its front side, position the switch on the left of the device body on position "O", being careful not to press the foot pedal. The device turns off.

PLEASE NOTE whenever the device is started, the following default setting is programmed: "function" ENDO, "irrigation" 3, "light" OFF.

05.2 → DESCRIPTION OF THE KEYBOARD

**Touch keyboard**

The user can configure the device by simply touching on the touch keyboard. Depending on the selected setting, the electronic feedback system automatically adjusts the correct operating frequency.

EN

**FUNCTIONS (ref. R inside cover)**

Depending on the type of surgery, it is possible to choose one of the 6 options available from the “function” list:

ENDO: dedicated to endodontic surgery and to the detachment of the sinus membrane


PERIO: dedicated to periodontal surgery

SPECIAL: dedicated only to the inserts for osteotomies with thicknesses of 0.35 mm and for prosthetic applications

CANCELLOUS: dedicated to the cutting and removal of poorly mineralized bone

CORTICAL: dedicated to the cutting and removal of highly mineralized bone

IMPLANT: dedicated to bone perforation in the technique of the implant site preparation

**FILLING AND CLEANING THE IRRIGATION CIRCUIT (ref. T inside cover)**

The device is equipped with the “pump/clean” key which, depending on the mode of use, allows to perform the PUMP function or the CLEAN function.

The PUMP function can be used at the beginning of the treatment, to flush the entire irrigation line up to the insert, so that the surgery can be started with the necessary irrigation (see paragraph 05.4).

The CLEAN function allows to perform a cleaning cycle of the irrigation circuit. This function must be performed when you have finished using the device and before cleaning, disinfecting, and sterilizing all the parts (see Cleaning and Sterilization Manual).

IRRIGATION (ref. Q inside cover)

The capacity of the peristaltic pump can be adjusted by selecting the numbers on the touch keyboard in the “irrigation” column.

7 capacity levels are foreseen:

0 = operation of the pump is closed: no irrigation outflows from the insert

From 1 to 6= the pump flow goes from 8 ml/min to approximately 75 ml/min.

The possibility of choosing the irrigation capacity levels is tied to the type of function selected, in the following manner:

ENDO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min)

PERIO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min)

SPECIAL - 6 flow levels: from 1 to 6



(from 8 to approximately 75ml/min)

CANCELLOUS - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)


CORTICAL - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)

IMPLANT - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)


PLEASE NOTE treatment without irrigation is possible only with the ENDO and PERIO functions, setting the irrigation capacity level on “0”.

In the event that an irrigation capacity lower than 8 ml/min is required, use the “bone grafting kit” (accessory that can be ordered separately), inserting it between the flow adjuster and the tube in silicon of the handpiece, making it pass through the peristaltic pump, and selecting 1 as the level of irrigation.


⚠ WARNING: If the “bone grafting kit” tube is kept inserted for the entire duration of the intervention, the capacity of the pump on all its levels is limited, independently of which insert is used.

LIGHT (ref. S inside cover)

Depending on the type of handpiece and operation that needs to be performed, it is possible to choose 3 options from the “light” list:

- By choosing the AUTO option, the LED light on the front terminal of the handpiece is switched on by pressing the foot pedal, and automatically switches off 3 seconds after the pedal is released.
- By choosing the ON option, the LED light on the front terminal of the handpiece stays on permanently, regardless of whether pressure is applied on the foot pedal or not. The light switches off 100 seconds after the last pressure of the foot pedal and the option shifts from ON to AUTO.
- By choosing the OFF option, the LED light on the front terminal of the handpiece stays off permanently.

The position of the LED light on the front terminal of the handpiece can be adjusted in the following way:

- Hold the body of the handpiece and lightly unscrew the metal ring nut located at the base of the front terminal, rotating it counter-clockwise.
- Rotate the front terminal so that the LED light goes into the desired and necessary position.
- To fasten it into position, screw the metal ring nut, rotating it clockwise.

SYMBOLS (ref. U inside cover)

The PIEZOSURGERY® *touch* is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol. To help the user identify the malfunctioning part, four symbols are foreseen which are described in paragraph 09.1.

05.3 → SAFETY PRECAUTIONS BEFORE AND DURING USE

⚠ WARNING: Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

⚠ CAUTION: Contraindication. Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ WARNING: Contraindications. Do not use the PIEZOSURGERY® *touch* on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ WARNING: Checking device status before the treatment. Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

⚠ CAUTION: CLEAN function. After the device is used with aggressive and non-aggressive solutions, it is necessary to perform a cleaning cycle on the tubes and the handpiece with the CLEAN function (see Cleaning and Sterilization Manual). If the tubes are not cleaned, the crystallization of the salts may seriously damage the device.

⚠ CAUTION: Contraindications. Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ WARNING: Infections control.
First use: All new and repaired accessories are supplied in NON STERILE conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilization Manual.
Subsequent uses: After every treatment, clean and sterilize all the reusable parts and accessories, following the instructions provided in the Cleaning and Sterilization Manual.

⚠ CAUTION: The electrical contacts inside the cord connector must be dry. Before connecting the handpiece to the device, make sure that the electrical contacts of the connector are perfectly dry, especially after the sterilization cycle in autoclave. If need be, dry the contacts by blowing compressed air onto them.

⚠ CAUTION: CLEAN function. The CLEAN function must be used after every treatment, before starting the cleaning and sterilization procedures.

⚠ WARNING: To ensure that the handpiece cools down, always activate it with the irrigation circuit correctly installed and filled. To fill the irrigation circuit, always use the PUMP function

⚠ WARNING: Treatments that require irrigation. Always check operation of the irrigation before and during use. Make sure the fluid outflows from the insert. Do not use the device if the irrigation does not work or if the pump is defective.

⚠ WARNING: Only use original Mectron inserts, accessories, and spare parts.

⚠ WARNING: Use of non-original Mectron inserts: this use entails finite damage to the handpiece threading, thus compromising correct operation and risking to cause harm to the patient.

ⓘ **CAUTION:** To correctly use the device, you must press the foot pedal and start it with the insert not making contact with the part to be treated, so that the electronic circuit can recognize the best point of resonance of the insert without interferences, allowing its optimal performance.

⚠ **WARNING:** Before every treatment, make sure that the insert appropriate for the treatment is inserted on the handpiece. Exclusively use the Mectron torque wrench to fasten the insert to the handpiece.

⚠ **WARNING:** The patient must not come into contact with the device body or the foot pedal.

⚠ **WARNING:** Do not change the insert while the handpiece is operating, to prevent causing wounds to the operator.

⚠ **WARNING: Breakage and wear-out of the inserts.** High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert. Do not bend, change shape of, or re-sharpen an insert in any way.

Bending an insert or applying leverage on it can lead to its breakage.

Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used. Excessive pressure on the inserts during their use can lead to their breakage. If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them.


The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent the patient from ingesting fragments of broken inserts.

When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out. Use of a worn-out insert reduces the cutting performance and can cause necrosis of the bone surface treated. During the intervention, frequently check that the insert is intact, especially its top.

During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use.

05.4 → INSTRUCTIONS FOR USE

After having connected all the accessories as illustrated in Paragraph 04.3, proceed as follows:

1


Open both the air socket on the flow adjuster and the liquid passage;

PIEZOSURGERY® touch

EN


To fill the irrigation circuit, use the PUMP function by selecting PUMP/CLEAN on the touch keyboard: all the other selection options present on the display are disabled, and the PUMP/CLEAN writing flashes;

2


While the PUMP/CLEAN wording is flashing, press the foot pedal once and release it; the PUMP/CLEAN stops flashing and the irrigation circuit starts to fill up;

3


As soon as the peristaltic pump starts, the entire scale of values of the "irrigation" section lights up and during the liquid passage, the value of the irrigation shifts from 6 to 0;

4


The cycle can be interrupted as soon you see liquid outflowing from the PIEZOSURGERY® touch handpiece, by pressing PUMP/CLEAN again or, alternatively, by pressing the foot pedal. The PUMP function is disabled and the keyboard is enabled again, and displays the last setting used;





5

Screw the chosen insert onto the PIEZOSURGERY® touch handpiece, until it reaches the limit stop;

6

7

Tighten the insert by using the Mectron torque wrench;

To correctly use the Mectron torque wrench, operate as follows:

8

Put the insert inside the wrench as shown;

9

Firmly hold the central body of the handpiece;

⚠ CAUTION: The handpiece must not be grabbed by its terminal part and/or cord, but only by its central body. The handpiece must not be rotated, but must be grasped firmly, and you must only rotate the wrench.

Rotate the wrench clockwise until the friction snaps (the external body of the wrench rotates compared to the handpiece, emitting mechanical “CLICK” sounds). The insert is now perfectly tightened;

10

On the keyboard, select the type of function and irrigation necessary and the light, if desired.

⚠ CAUTION: For the correct setting of the Function parameters with the specific insert you wish to use, consult the Chart annexed to this manual titled “Appropriate settings for the inserts on the PIEZOSURGERY® touch” or the illustrative leaflet of the Mectron insert you’ve purchased.

05.5 → IMPORTANT INFORMATION ON THE INSERTS

WARNING:

- When the titanium nitride coating is visibly worn out, the insert must be replaced. Use of an overly worn out insert reduces its cutting efficiency.
- Diamond-coated inserts: the diamond inserts must be replaced when the titanium nitride coating is visibly worn out and, at any rate, after a maximum of 10 treatments.
- Do not activate the handpiece while the insert is in contact with the treated part, so that the electronic circuit can recognize the best point of resonance of the insert and allow its optimal performance.
- Check the condition of wear of the insert and that it is intact before and during every use. In the event that a drop in performance becomes apparent, proceed to replace it.
- Use original Mectron inserts only. Use of non-original inserts, in addition to voiding the warranty, damages the threading of the PIEZOSURGERY® *touch* handpiece, with the risk of no longer being able to screw the original inserts correctly during subsequent use. Moreover, the device settings are tested and guaranteed to operate correctly only when original Mectron inserts are used.
- Do not change the shape of the insert in any way by either bending or filing it. This could cause it to break.
- Do not use an insert that has suffered any type of deformation.
- Do not attempt to sharpen the insert used.
- Always check that the threaded parts of the insert and of the handpiece are perfectly clean – see the Cleaning and Sterilization Manual.
- If excessive pressure is applied on the insert, it can cause the insert to break and possibly harm the patient.
- For information on how to correctly use the inserts, consult the annexed sheet “Appropriate settings for the inserts on the PIEZOSURGERY® *touch*” or the illustrative leaflet of the Mectron insert you’ve purchased.
- Before using the PIEZOSURGERY® *touch*, make sure you have prepared the operatory site by having first moved away the soft tissues, to avoid damaging them. It may happen that, while cutting the bone, accidental contact of certain parts of the insert with the soft tissues inflicts small traumas. To minimize this risk, use specific protective instruments.

06 → MAINTENANCE

If the device is not used for prolonged time, observe the following recommendations:

- 1 Run a complete cleaning cycle on the irrigation circuit with the CLEAN function (see the Cleaning and Sterilization Manual);
- 2 Disconnect the device from the electrical network;
- 3 If the period of disuse is prolonged, put the device back in its original package and store it in a safe place;
- 4 Prior to using the device again, clean and sterilize the handpiece, the inserts, the

wrench, the tubes, and the connectors, following the instructions provided in the Cleaning and Sterilization Manual;

- 5 Check that the inserts are not worn out, deformed, or broken, placing special attention to the integrity of their tip.

⚠ WARNING: Periodically check that the electrical power cable is intact; if it is damaged, replace it with an original Mectron spare part.

EN

07 → DISPOSAL MODES AND PRECAUTIONS

⚠ WARNING: Hospital wastes.

Hospital wastes. Treat the following items as hospital waste:

- Inserts: when they are worn out or broken;
- Flow adjuster: at the end of each intervention;
- Peristaltic pump tube: after 8 sterilization cycles;
- Torque wrench for inserts: when worn out or broken.

The use and throw-away materials and materials that entail a biological risk must be disposed of in accordance to the norms in force locally concerning hospital wastes.

PIEZOSURGERY® *touch* must be disposed of and treated as a waste for separate collection. Disregard of the previous points may entail a fine, pursuant to Directive 2002/96/CE. It is up to the purchaser to hand over the device for its disposal to the retailer who supplies him new equipment; the instructions for proper disposal are available from Mectron.

08 → TECHNICAL DATA

EN

Device compliant to Dir. 93/42/CEE:	Class II a
Classification as per EN 60601-1:	I Applied part type B (handpiece, insert) IP 20 (device) IP X8 (foot pedal)
Device for intermittent operation:	60sec. ON - 30sec. OFF with irrigation 30sec. ON - 120sec. OFF without irrigation(ENDO,PERIO)
Power supply voltage:	100-240 Vac 50/60 Hz
Max. power absorbed:	120 VA
Fuses:	Type 5 x 20 mm T 2AL, 250V
Operating frequency:	Automatic scan From 24 KHz to 36 KHz
Power types:	ENDO PERIO SPECIAL CANCELLOUS CORTICAL IMPLANT
Peristaltic pump capacity:	Adjustable on the touch screen: ENDO / PERIO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min) SPECIAL / CANCELLOUS / CORTICAL / IMPLANT 6 capacity levels: from 1 to 6 (from 8 to approximately 75ml/min)
LED system of the handpiece:	Light function on AUTO: LED system of the handpiece switches on as soon as the device starts to operate and switches off 3 seconds after the foot pedal is released Light function on ON: The LED of the handpiece is always on; after 100 seconds of foot pedal disuse it switches off by itself and the light function moves into position AUTO Light Function on OFF: The LED of the handpiece is always off.
Protections of the APC circuit:	No handpiece detected Cord interruption Insert not tightened correctly or broken
Operating conditions:	from +10°C to +40°C Relative humidity from 30% to 75% Pressure of air P: 800hPa/1060hPa
Transport and storage conditions:	from -10°C to +70°C Relative humidity from 10% to 90% Pressure of air P: 500hPa/1060hPa
Peristaltic pump tube:	It is advisable not to exceed 8 sterilization cycles
Weights and sizes:	3,2 Kg L -W - H 300 x 250 x 95 mm

08.1 → ELECTROMAGNETIC COMPATIBILITY EN 60601-1-2

⚠ WARNING: Interference with other equipment. Though compliant with the standard IEC 60601-1-2, the PIEZOSURGERY® *touch* may nonetheless interfere with other devices nearby. The PIEZOSURGERY® *touch* must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

⚠ WARNING: Portable and mobile radio communication appliances may affect the correct functioning of the device.

⚠ WARNING: Interference from other equipment. An electrical scalpel or other electro-surgical units near the PIEZOSURGERY® *touch* device may interfere with its correct operation.

⚠ WARNING: The device requires specific EMC precautions and must be installed and activated in accordance with the EMC information given in this paragraph.

Guidance and manufacturer's declaration - Electromagnetic emissions

The PIEZOSURGERY® *touch* is intended for use in the electromagnetic environment specified below. The customer or user of the PIEZOSURGERY® *touch* should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The PIEZOSURGERY® <i>touch</i> only uses RF energy for internal function. Therefore, its FR emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The PIEZOSURGERY® <i>touch</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - Electromagnetic immunity

The PIEZOSURGERY® touch is intended for use in the electromagnetic environment specified below. The customer or user of the PIEZOSURGERY® touch should ensure that it is used in such an environment.

EN

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	The device continues to work regularly and in safety	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	The device continues to work regularly and in safety	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	The device continues to work regularly and in safety	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycles 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	The device can vary from the required levels of immunity with a duration of <5% / >95% / 5 s as long as the device remains in safety, no malfunctions have been detected and can be restored to pre-test status with the intervention of the operator	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	The device continues to work regularly and in safety	Power frequency magnetic fields should be at levels characteristic of a typical location in a commercial or hospital environment.

N.B.: U_T is the AC mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - Electromagnetic immunity

The PIEZOSURGERY® *touch* is intended for use in the electromagnetic environment specified below.

The customer or user of the PIEZOSURGERY® *touch* should ensure that it is used in such an environment.

EN

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment Guidance
Conducted RF IEC 61000-4-6	3 Veff from 150 kHz to 80 MHz	The device continues to work regularly and in safety	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m from 80 MHz to 2,5 GHz		Recommended separation distance $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ from 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ from 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , may be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

N.B.:

(1) at 80 MHz and 800 MHz, the higher frequency range applies.

(2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with any reasonable accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PIEZOSURGERY® *touch* is used exceeds the applicable RF compliance level given above, the PIEZOSURGERY® *touch* should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the PIEZOSURGERY® *touch*.

b Over the frequency range from 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the PIEZOSURGERY® touch

The PIEZOSURGERY® touch is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the PIEZOSURGERY® touch can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile FR communications equipment (transmitters) and the PIEZOSURGERY® touch as recommended below, according to the maximum output power of the communications equipment.

EN

Rated maximum output power of transmitter 'W'	Separation distance according to the frequency of transmitter 'm'		
	from 150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	from 80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	from 800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be calculated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

N.B.:

- (1) at 80 MHz and 800 MHz, the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

09 → TROUBLESHOOTING

09.1 → DIAGNOSTIC SYSTEM AND SYMBOLS ON KEYBOARD

II PIEZOSURGERY® *touch* is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol. By using the following chart, the user is guided toward the identification and possible solution of the malfunction detected.

EN

Symbols on keyboard	Possible cause	Solution
	Handpiece cord contacts wet	Thoroughly dry the contacts with compressed air
	PIEZOSURGERY® <i>touch</i> handpiece not connected to device	Connect the handpiece
	Handpiece failure	Replace the handpiece
	Sync circuit malfunction	Contact an Authorized Mectron Service Center
	Insert not tightened correctly on handpiece	Unscrew the insert and correctly screw it on again with the Mectron torque wrench (see paragraph 05.4)
	Insert broken, worn-out or deformed	Replace the insert
	Handpiece cord contacts wet	Thoroughly dry the contacts with compressed air
	Peristaltic pump malfunction	Check that there are no impediments to pump rotation
	Silicon tubing not positioned correctly inside the pump	Correctly reposition the silicon tubing inside the pump (see paragraph 04.3)
	The device has been turned off and on again without waiting 5 seconds	Turn device off and wait 5 seconds before turning it on again
	Abnormalities on electrical network or excessive electrostatic discharges or internal abnormalities	Turn device off and wait 5 seconds before turning it on again If the signal persists, contact an Authorized Mectron Service Center
	Turn-on procedure incorrect: the device has been turned on with the foot pedal pressed	Check that the foot pedal is not pressed. If the problem persists, disconnect the pedal and, if need be, contact an Authorized Mectron Service Center

09.2 → QUICK SOLUTION TO PROBLEMS

EN

Problem	Possible cause	Solution
The device does not turn on after having brought the switch into position “1”	The electrical power cable terminal is badly inserted in the rear plug of the device	Check that the power supply cable is firmly connected
	The electrical power cable is defective	Check that the power supply socket works properly. Replace the electrical power cable
	The fuses are out of order	Replace the fuses (see paragraph 09.3)
The device is on but not working. The display does not signal any error	The foot pedal plug is incorrectly inserted in the device socket	Correctly insert the pedal plug in the socket on the back of the device
	The foot pedal does not work	Contact an Authorized Mectron Service Center
The device is on but not working. One of the following symbols appears on the screen: 	See paragraph 09.1 for the possible cause, according to the symbol that has been displayed	See paragraph 09.1 for the action to undertake, according to the symbol that has been displayed
A slight whistling sound coming from the PIEZOSURGERY® touch handpiece is heard during operation.	The insert is not correctly tightened on the handpiece	Unscrew and correctly screw the insert again with the Mectron torque wrench (see paragraph 05.4)
	The irrigation circuit has not been completely filled	Fill the irrigation circuit via the PUMP function (see paragraph 05.4)

Problem	Possible cause	Solution
No liquid outflows from the insert during operation	The insert type does not foresee liquid passage	Use an insert type with liquid passage
	The insert is obstructed	Unscrew the insert from the handpiece and free the insert water passage by blowing compressed through it. If the problem persists, replace the insert with a new one
	The handpiece is obstructed	Contact an Authorized Mectron Service Center
	The irrigation level on the screen is adjusted on "0"	Adjust the irrigation level
	The liquid bag is empty	Replace the bag with a full one
	The air inlet of the flow adjuster has not been opened	Open the air inlet of the flow adjuster
	The silicon tubing are incorrectly installed	Check the connections of silicon tubing
The device works properly but the pump strains/ stresses	Excessive pressure of the impeller on the peristaltic pump tube	Check that the tube inside the peristaltic pump is properly positioned (See paragraph 04.3)
The pump turns correctly, but when it stops, liquid outflows from the handpiece	The lid of the peristaltic pump is not closed correctly	Check that the lid of the peristaltic pump is perfectly closed (See paragraph 04.3)
Insufficient performance	The insert is not correctly tightened on the handpiece	Unscrew and correctly screw the insert again with the Mectron torque wrench
	Insert broken, worn-out, or deformed	Replace the insert with a new one


09.3 → REPLACEMENT OF THE FUSES

⚠ WARNING: Switch the device off.
Always turn the device off with the main switch and disconnect it from the electrical power socket before performing the following intervention.

EN

Apply leverage with a flat screwdriver, inserting its tip in the seat of the fuse-holder drawer located under the power supply socket;


1

Pull out the fuse-holder drawer;


2

⚠ WARNING: Replace the fuses, complying to the characteristics indicated in Chapter 08 - TECHNICAL DATA

Reinsert the drawer in its housing.

3

09.4 → SENDING THE DEVICE TO AN AUTHORISED MECTRON SERVICE CENTRE

In the event that technical assistance is needed on the device, contact one of the Authorized Mectron Service Centers or your Retailer. Do not attempt to repair or modify the device and its accessories.

Clean and sterilize all the parts that need to be sent to an Authorized Mectron Service Center, following the instructions provided in the Cleaning and Sterilization Manual supplied with the device.

Leave the sterilized parts in the packet, procedure that confirms the sterilization process has been performed.

The cleaning and sterilization requirements are in line with those in force concerning the safeguarding of health and safety on the workplace, as per Law Decree 626/94 and DLgs 81/08 and its subsequent amendments, both laws of the Italian State.

In the case that the customer does not fulfill the requirements, Mectron reserves the right to charge him or her the cleaning and sterilization expenses, or to reject the goods received in unsuitable conditions, returning them to the customer, at his or her expense, for them to be correctly cleaned and sterilized.

The device must be returned properly packaged, accompanied by all the accessories and by a sheet bearing the:

- Data of the owner with telephone number
- Product name
- Serial number and/or lot number
- Reason for goods returned / description of the malfunction
- Photocopy of delivery note or purchase invoice of the device

CAUTION: Package

Pack the device in its original package to prevent damages during transport.

Once the material has been received by the Authorized Mectron Service Center, qualified technical personnel will provide its evaluation to the given circumstance. The repair works will be performed only upon their prior acceptance by the end customer. For further details, contact the closest Authorized Mectron Service Center or your retailer.

Unauthorized repair works may damage the system and make the warranty void, and absolve Mectron of any responsibility for direct or indirect harm to persons or damages to things.

10 → WARRANTY

EN

All the Mectron devices are subjected to an accurate end inspection that ascertains their full functionality before they are placed on the market.

Mectron provides warranty on the PIEZOSURGERY® touch purchased new from a retailer or a Mectron importer to cover material or manufacturing defects for:

- 2 YEARS (TWO) on the device from the date of purchase;
- 1 YEAR (ONE) on the handpiece complete with cord from the date of purchase.

The accessories are not included in the warranty.

During the period in which the warranty is valid, Mectron commits itself to repair (or at its own free choice, replace) those product parts that were to prove defective according to its judgment, free of charge.

Entire replacement of the Mectron products is excluded.

The warranty of the manufacturer and the homologation of the device are not valid in the following cases:

- The device is not used in accordance to the intended use foreseen for it.
- The device is not used in compliance to all the instructions and precautions described in this manual.
- The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions.
- The assembly operations, extensions, adjustments, updates, and repair works are performed by personnel not authorized by Mectron.

- The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 08 - TECHNICAL DATA.
- Use of non-original Mectron inserts, accessories, and spare parts that can compromise correct device operation and cause harm to the patient.
- Accidental breakage due to transport
- Damages due to incorrect use or negligence, or to connection to a power voltage other than the one foreseen
- warranty expiry

PLEASE NOTE The warranty is valid only if the warranty sheet annexed to the product has been filled out in all its parts, and sent to our offices or possibly to the retailer or Mectron importer within 20 (twenty) days from the date of purchase, which is the date that appears on the delivery note/purchase invoice issued by the retailer/importer.

To benefit of the warranty service, the customer must return, at his own expense, the device to be repaired to the retailer/Mectron importer from whom he has purchased the product.

See paragraph 09.4 for the details relative to the procedure to send the device to an Authorized Mectron Service Center.

The indications that appear in this publication are not binding and can be modified without fore-notice.

The Italian version of this manual is the original document from which its translations have been obtained. In case of any discrepancy, the Italian version will have pertinence.

Texts, images, and graphics of this manual are property of Mectron S.p.A., Carasco, Italy. All rights reserved.

The contents cannot be copied, distributed, changed, or made available to third parties without the written approval of Mectron S.p.A.

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→ SINUS LIFT TECHNIQUE

Sentineri R., Dagnino G.

Sinus Augmentation by crestal approach with the Sinus Physiolift device
Implants. 2011; 1:13-17

Vercellotti T, De Paoli S, Nevins M.

The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure

Int J Periodontics Restorative Dent. 2001; 21(6):561-567

Vercellotti T, Nevins M, Jensen Ole T.

Piezoelectric Bone Surgery for Sinus Bone Grafting.

The Sinus Bone Graft, Second Edition. Edited by Ole T. Jensen, Quintessence Books. 2006; 23:273-279

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Wallace SS, Mazor Z, Froum SJ, Cho SC, Tarnow DP.

Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.

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Muñoz-Guerra MF, Naval-Gías L, Capote-Moreno A.

Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®.

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Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May; 107(5):648-55.

→ RIDGE EXPANSION

Vercellotti T.

Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique.

Int J Periodontics Restorative Dent. 2000; 20(4):359-365

Vercellotti T, Russo C, Gianotti S

A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.

World Dentistry. 2000 online article

Enislidis G, Wittwer G, Ewers R.

Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.

Int J Oral Maxillofac Implants. 2006; 21(3):445-449

→ CORTICOTOMY TECHNIQUE

Grenga V., Bovi M.

Corticotomy-Enhanced Intrusion of an Overerupted Molar Using Skeletal Anchorage and Ultrasonic Surgery.

JCO, Inc..2013; Januara:50-55

Vercellotti T, Podestà A.

Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement.

Int J Periodontics Restorative Dent. 2007; 27:325-331

→ BONE GRAFTING

Chiriac G, Herten M, Schwarz F, Rothamel D, Becker J.

Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation.

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Implant. 2005; 11(4):261-274

Sivolella S, Berengo M, Scarin M, Mella F, Martinelli F.

Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study.

Arch Oral Biol. 2006; 51(10):883-891

Happe A.

Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases.

Int J Periodontics Restorative Dent. 2007; 27(3):241-249

EN

IT

DE

FR

ES

SV

Patel A, Schofield J.

Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation.

Implant Dentistry Today. 2007; 1(4):20-24

Sohn DS, Ahn MR, Lee WH, Yeo DS, Lim SY.

Piezoelectric osteotomy for intraoral harvesting of bone blocks.

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Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report.

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Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study.

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Schweiz Monatsschr Zahnmed. 2005;115(8):692-703.

Gellrich NC, Held U, Schoen R, Pailing T, Schramm A, Bormann KH.

Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures.

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A comparative study of bone densitometry during osseointegration: Piezoelectric surgery versus rotatory protocols.

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Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.

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Changes in Implant Stability Using Different Site Preparation Techniques: Twist Drills versus Piezosurgery. A Single-Blinded, Randomized, Controlled Clinica Trial.

Clinical Implant Dentistry and Related Research, 2011:1-10

Bengazi F., Lang N.P., Canciani E., Viganò P., Urbizo Velez J., Botticelli D.

Osseointegration of implants with dendrimers surface characteristics installed conventionally or with Piezosurgery®. A comparative study in the dog.

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Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.

J Periodontol. 2007; 78(4):716-722

→ OSTEOTOMY CLOSE TO NERVES

Bovi M.

Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report.

Int J Periodontics Restorative Dent. 2005; 25(4):375-383

Geha H, Gleizal A, Nimeskern N, Beziat JL.

Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery.

Plast Reconstr Surg. 2006; 118(7):1598-1607

Schaeren S, Jaquière C, Heberer M, Tolnay M, Vercellotti T, Martin I.

Assessment of nerve damage using a novel ultrasonic device for bone cutting.

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J Oral Maxillofac Surg. 2008; 66(3):593-596

Sakkas N, Otten JE, Gutwald R, Schmelzeisen R.

Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report.

Br J Oral Maxillofac Surg. 2008; 46(4):270-271

Tordjman S, Botoli LT.

Implants Juxta- canauxaires

L'information Dentaire, 2007 May;89(26):1499.

Degerliyurt K, Akar V, Denizci S, Yucel E.

Bone lid technique with PIEZOSURGERY® to preserve inferior alveolar nerve.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 Dec;108(6):e1-5.

→ EXTRACTIONS

Grega V, M. Bovi.

Piezoelectric Surgery for Exposure of Palatally Impacted Canines.

J Clin Orthod. 2004; 38(8):446-448

Sivolella S, Berengo M, Fiorot M, Mazzuchin M.

Retrieval of blade implants with piezosurgery: two clinical cases.

Minerva Stomatol. 2007; 56(1-2):53-61

EN

IT

DE

FR

ES

SV

Sortino F, Pedullà E, Masoli V.

The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery.

J Oral Maxillofac Surg. 2008 Dec;66(12):2444-8

→ **PERIODONTAL SURGERY**

Vercellotti T, Pollack AS.

A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery.

Compend Contin Educ Dent. 2006; 27(5):319-325

Seshan H, Konuganti K, Zope S

PIEZOSURGERY® in periodontology and oral implantology.

J Indian Soc Periodontol. 2009 Sep;13(3):155-6

→ **DISTRACTION OSTEOGENESIS**

González-García A, Diniz-Freitas M, Somoza-Martin M, García-García A.

Piezoelectric Bone Surgery Applied in Alveolar Distraction Osteogenesis: A Technical Note.

Int J Oral Maxillofac Implants. 2007; 22(6):1012-1016

Menini I, Zornitta C, Menini G.

Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report.

Int J Periodontics Restorative Dent. 2008; 28(2):189-196

Lee HJ, Ahn MR, Sohn DS.

Piezoelectric distraction osteogenesis in the atrophic maxillary anterior area: a case report.

Implant Dent. 2007 Sep;16(3):227-34.

→ **HISTOLOGICAL ADVANTAGES**

Vercellotti T, Nevins ML, Kim DM, Nevins M, Wada K, Schenk RK, Fiorellini JP.

Osseous Response following Resective Therapy with a Piezosurgery®.

Int J Periodontics Restorative Dent. 2005; 25(6):543-549

Gleizal A, Li S, Pialat JB, Béziat JL.

Transcriptional expression of calvarial bone after treatment with low - intensity ultrasound: An in vitro study.

Ultrasound Med Biol. 2006; 32(10):1569-1574

Stübinger S, Goethe JW.

Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications

Journal of Oral and Maxillofacial Surgery 2007 Sep;65(9):39.e7-39.e8.

Sohn DS, Lee JK, An KM, Shin HI.

Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report.

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Inf Dent. 2004; 86(41):2887-2893

Vercellotti T.

Technological characteristics and clinical indications of piezoelectric bone surgery.

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Schweiz Monatsschr Zahnmed. 2004; 114(4):365-377

Troiani C, Russo C, Ballarani G, Vercellotti T.

Piezoelectric Surgery: A new reality to cut and manage bone.

Maxillo Odontostomatologia- International Journal of Maxillo Odontostomatology - S.I.M.O. 2005; 4(1):23-28

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Z Zahnärztl Impl. 2005; 21(1):48-59

Schlee M, Steigmann M, Bratu E, Garg AK.

Piezosurgery: basics and possibilities.

Implant Dent. 2006; 15(4):334-340

Sembronio S, Albiero AM, Polini F, Robiony M, Politi M.

Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007; 104(1):e7-e10

Su YC.

Development and clinical application of ultrasonic osteotomy in dentistry.

Shanghai kou qiang yi xue = Shanghai journal of stomatology. 2007; 16(1):1-7 [chinese]

Peivandi A, Bugnet R, Debize E, Gleizal A, Dohan DM.

Piezoelectric osteotomy: applications in periodontal and implant surgery.

Rev Stomatol Chir Maxillofac. 2007; 108(5):431-440

Labanca M, Azzola F, Vinci R, Rodella LF

Piezoelectric surgery: Twenty years of use.

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Revue Implantologie, 2006 Feb;5-18.

➔ OVERVIEW ON MAXILLOFACIAL SURGERY APPLICATIONS

Robiony M, Polini F, Costa F, Vercellotti T, Politi M.
Piezoelectric Bone Cutting in multipiece
maxillary osteotomies. Technical Note.

J Oral Maxillofac Surg. 2004; 62:759-761.

Eggers G, Klein J, Blank J, Hassfeld S.
PIEZOSURGERY®: an ultrasound device for
cutting bone and its use and limitations in

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Beziat JL, Béra,JC, Lavandier B, Gleizal A.
Ultrasonic osteotomy as a new technique in
craniomaxillofacial surgery.

International Journal of Maxillo-facial Surgery,
2007;36(6):493-500.

Gleizal A, Béra JC, Lavandier B., Béziat JL.
Piezoelectric osteotomy: a new technique
for bone surgery – advantages in craniofacial
surgery.

Childs Nerv Syst. 2007;23(5):509-513.

Guo ZZ, Liu X, Li Y, Deng YF, Wang Y
The use of PIEZOSURGERY® osteotomy
in treatment of longstanding maxillary
fractures: report of 12 consecutive patients.

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G Chir 2007 May; 28(5):222-6.

Robiony M, Polini F, Costa F, Zerman N, Politi M.
Ultrasound bone cutting for surgically
assisted rapid maxillary expansion (SARME)
under local anaesthesia.

Int J Oral Maxillofac Surg. 2007;36(3):267-9.

Landes CA, Stübinger S, Rieger J, Williger B, Ha TK, Sader R.
Critical evaluation of piezoelectric osteotomy
in orthognathic surgery: operative technique,
blood loss, time requirement, nerve and
vessel integrity.

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Rev Stomatol Chir Maxillofac. 2008 Dec;109(6):375-8.

Landes CA, Stübinger S, Ballon A, Sader R
Piezoosteotomy in orthognathic surgery
versus conventional saw and chisel
osteotomy.

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