

RS wash

Perfect finish for dental products from the printer: automatic and environmentally-friendly cleaning.

The RS wash automated cleaning system boasts several new features to make cleaning particularly safe and easy (patents pending). Fully automatic cleaning in just about 6-8 minutes and clean handling of sticky resins are made possible. The consumption of cleaning agent is reduced to a minimum, while ensuring optimal results.

For easiest handling, the RS wash is equipped with a data link to the printer to preselect the correct cleaning program and the matching cleaning agent.

Performance parameters	RS wash
Volume	$130 \times 75 \times 60 \text{ mm}$
Cleaning time	approx 6–8 minutes (depending on material)
Cleaning medium	Isopropanol, metylalcohol, anorganic, (depending on material)
Connections	WLAN/LAN
Dimension (W \times H \times D)	230 imes 270 imes 450 mm

RS wash	
	r
rapidshape	odd
_	ropidshop





sticky resins



Intelligent connectivity to printer

No time consuming cleaning