








A YOUNGER SMILE, NOW.

Professional Teeth Whitening

+ DESENSITIZING KIT



syringe of LWS desensitizing gel



BEFORE



AFTER

FAST AND SURPRISING RESULTS.

The laser is the fastest and most effective way of achieving naturally white teeth, without causing any damage. The oxidization process set off by the hydrogen peroxide inside the **Doctor Smile LWS gel** is activated by laser light after being spread on the surface of teeth. As it dissolves it releases free radicals, and thanks to the interaction of the laser and the peroxide with the molecular pigments, creates a whitening effect. Doctor Smile LWS TITANIUM contains the pigment TiO₂ (titanium dioxide), a photocatalytic molecule that is naturally white, and filters the light to avoid any overheating. The presence of titanium dioxide allows the use of a minor concentration of hydrogen peroxide in order to maximize the bleaching effect and minimize side effects.



PRODUCT CODE L P LWS 003.1

DISTRIBUTED BY DOCTOR SMILE



MORE INFO

doctor-smile.com
info@doctor-smile.com

The LWS Laser Whitening **System** is a professional bleaching product used in combination with laser wavelength (diode) to 1064nm (Nd:YAG).

DOCTOR SMILE IS A TRADEMARK OF LAMBDA SPA

Via dell'Impresa, 1 36040 BRENDOLA (VI) Italy T +39 0444.349165 F+39 0444.349954