

VRN

DS-100

Ultrasurgery Piezo Bone Surgery

Efficient, safe and minimally invasive



Features

7 inch large LED touch screen

Doctor-friendly UI

Multiple working modes

LED surgery handpiece

Intelligent recognition of soft tissue


Separate liquid path design

Multi-function foot pedal

<https://stomshop.pro>

DS-100

Ultrasurgery Piezo Bone Surgery



Technical parameter

Power supply voltage:	100-240V 50/60Hz 120VA
Device for intermittent operation:	60s ON, 10s OFF
Fuses:	T2AH 250V
Flow:	25-100ml/min
Radial amplitude:	8-63um
Tranverse amplitude:	1.5-12um
Working tips frequency:	22.3kHz-29.6kHz
Waterproof level:	IPX0 (device), IPX1 (foot pedal)
Tips (basic kit):	US1, US2, US4, US5, UL3, UC1
Size of main unit:	321×310×138 (L×W×H)
Weight of main unit:	3.8KG