

# Ez SHOT HD

Easy to Use

Auto Focus HD intraoral camera



Wired type

-IC-WHT60 -IC-WHW61  
 -IC-DCAM70 -IC-DCAU71  
 -IC-WHC50 -IC-EZAM80  
 -IC-EZAU81

-CL-DLT30 -CL-DL10  
 -CL-DC20 -CL-AT24

-AL-DF10

# Easy to Use

## Auto Focus HD intraoral camera

### 1 Features of Ez SHOT HD

#### Wireless type only

- High speed digital wireless**  
High speed digital wireless provide more clear live image without interruption.
- Wireless receiver**  
Compact wireless receiver without any power supply.

- Macro to Full arch**  
Camera is able to focus from macro mode (close) to full arch.
- Auto focus**  
Can get clear image with auto focus function in oral cavity.
- HD image**  
Can get high definition image without distortion.
- Shadowless Lighting**  
Rounded 12pcs LEDs make shadowless lighting.
- Adjustable LED brightness**  
To avoid image quality loss by light reflection, can adjust LED brightness.
- Black & White function**  
Can take and send analog X-ray film image to PC.
- Plug & Play**  
No need additional driver to install the device.
- Easy to capture**  
Easy to image capture by camera only without extra device like foot pedal.
- Software compatibility**  
Can be use camera with most of dental imaging software.  
\* It would need twain driver for some of dental imaging software.




### 2 Specification

Model name	IC-EZAM80(Wireless model)	
Camera	Operation voltage	3.7V
	Size	Approx. 255 x 27 x 29 mm
	Weight	Approx. 89g
	Resolution	HD(720P)
	Sensor	1.3M CMOS
Charger	Operation voltage	5V
	Size	Approx. 110 x 105 x 68 mm
	Weight	Approx. 100g
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

Model name	IC-EZAU81(Wired model)	
Camera	Operation voltage	5V
	Size	Approx. 202 x 27 x 29 mm
	Weight	Approx. 55g
	Resolution	HD(720P)
	Sensor	1.3M CMOS
USB cable	Length	Approx. 3M
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

### 3 Minimum System Requirement

- PC system: 2.0 Ghz CPU, 1GB RAM
- OS: Windows XP service pack 3
- USB port: USB 2.0



Upper images are take by real product.

### 4 DrsViewer2

Image Viewer included.

