

# OPTIMIZED FOR YOUR DAILY APPLICATIONS

## Identica *light*

The most **COST EFFECTIVE**, **RELIABLE**, and **CONVENIENT WAY**  
for all your dental applications



# PROVEN RELIABILITY. Identica *light*

## Increase your work productivity with Identica Light!

Identica Light is optimized for your daily dental applications, guaranteeing reliable scan results with the minimum investment of time and money. You can now digitize your dental work-flow and get fast and improved scanning ability but at a very affordable price.

Furthermore, Identica light is able to export the scan data in the open STL format that works with any CAD software.

### Key features



#### OPEN SPACIOUS

The Identica Light has a fully open platform where conventional closed-platform scanners do not, bringing the maximum in user convenience to your scanning solutions for all your dental applications.



#### COST EFFECTIVE

You can achieve a digitized dental laboratory with the affordable and reasonably priced Identica Light.



#### INTUITIVE SCAN PROCESS

The Identica Light comes with scanning software that employs very simple and easy to follow scan sequences.



#### OPEN SYSTEM

The basic output format of the Identica Light is STL, so any CAD software can utilize the file.



### Specification

Category	Description
Light Source/Color	LED, 50 ANSI-lumens / Cyan
Scanning Area	80mm x 60mm x 60mm
Dimension	290mm x 290mm x 340mm
Weight	16 kg
Interface	USB 3.0 B Type
Power	AC 100~240V, 50~60 Hz