

COMPATIBLE MATERIAL LIST

SEPTEMBER 2019

V E R S V S
— 3 8 5 —

DENTAL 3D PRINTER



DETAX

Ceramics with passion

	50µm layer height	100µm layer height	Application
Freeprint® Model Grey	✓	✓	Model
Freeprint® Model T	✓	✓	Model for thermoforming
Freeprint® Temp A1	✓	✓	Temporary crown and bridge
Freeprint® Temp A2	✓	✓	Temporary crown and bridge
Freeprint® Temp A3	✓	✓	Temporary crown and bridge
Freeprint® Ortho UV	✓	✓	Orthodontic appliances
Freeprint® Cast 2.0	✓	✓	Metal casting
Freeprint® Cast	✓	✓	Metal casting
Freeprint® Model Ivory	✓	✓	Model
Freeprint® Model Sand	✓	✓	Model
Freeprint® Denture	✓	✓	Denture base
Freeprint® Gingiva	✓	✓	Gingiva Mask
Freeprint® Tray	✓	✓	Individual tray



OptiPrint Model GoldenBrown	✓	✓	Model
OptiPrint Model Ivory (385nm)	✓	✓	Model
OptiPrint Guide (385nm)	✓	✓	Splint, Surgical Guide
OptiPrint Cast	✓	✓	Metal casting
OptiPrint Gingiva	✓	✓	Gingiva mask
OptiPrint Temp A2 385	✓	✓	Temporary crown and bridge
OptiPrint Temp A3 385	✓	✓	Temporary crown and bridge
OptiPrint Denture Dark Pink 385	✓	✓	Denture Base
OptiPrint Denture Light Pink 385	✓	✓	Denture Base
OptiPrint Temp B2 385	✓	✓	Temporary crown and bridge
OptiPrint MockUp	✓	✓	Try-In
OptiPrint Match	✓	✓	Model scan
OptiPrint Tray Neon Green	✓	✓	Individual tray



C&B A3.5	✓	✓	Temporary crown and bridge
SG	✓	✓	Surgical guide
Model 2.0	✓	✓	Model
Base	✓	✓	Denture base
Gingiva	✓	✓	Gingiva Mask
Try-In	✓	✓	Try-ins
Ortho Rigid	✓	✓	Splints
Ortho Clear	✓	✓	Splints



VeriModel OS White	✓	✓	Model
VeriModel OS Ivory	✓	✓	Model
VeriGuide OS Clear	✓	✓	Splints, Surgical guide
VeryTray Neon	✓	✓	Custom tray
VeriGuide	✓	✓	Splints, Surgical guide



S-Dent GT Medical	✓	✓	Model
S-Dent Ocre MTD	✓	✓	Model
S-Cast	✓	✓	Metal casting



FotoDent® Model	✓	✓	Model
FotoDent® Guide 385 nm	✓	✓	Surgical guide
FotoDent® Tray			Cutoms tray
FotoDent® Cast	✓	✓	Metal casting
FotoDent® Gingiva 385 nm	✓	✓	Gingiva mask
FotoDent® Model2 385 nm	✓	✓	Model

FUNTODO

Dentifix	✓	✓	Model
Snow White	✓	✓	Model



KeyModel	✓	✓	Model
KeyOrthoModel	✓	✓	Model
KeyMask	✓	✓	Gingiva mask
KeyGuide	✓	✓	Guides
KeySplint Hard	✓	✓	Splints
KeySplint Soft	✓	✓	Soft splints
KeyCast	✓	✓	Metal casting



GC Temp Print (Composite) Light	✓		Class IIA Composite for temporary
GC Temp Print (Composite) Medium	✓		Class IIA Composite for temporary

DeltaMed

TURNING IDEAS INTO MATERIALS

3Delta Model 320	✓	✓	Model
3Delta Model Ortho	✓	✓	Fast Ortho Model
3Delta Model Silbergrau	✓	✓	Model
3Delta Cast P	✓	✓	Ceramic pressing casting
3Delta Cast M	✓	✓	Metal casting