

## With Swinging instruments hoses :

It's a new generation unit, a combination of the high Swiss technology with Italian production and design. Ergonomic handling, smooth operating and a high quality at very competitive price.



## The Electromechanical Patient Chair "AG04SW" with motor lifting to low tension

- Doted with trendelenburg synchronized movement
- The headrest is provided with a double articulation
- 4=3 different chair programs positions and plus one additional rising and return to the last work position
- All positions are stored and controlled from the dentist side and part from the assistant's side
- The chair is equipped with 2 safety switches
- No cables or supply hoses on floor

## The Hydro Unit :

- Removable and auto-cleavable ceramic cuspidor bowl
- Automatic timer of cup filler and bowl flush.
- Tube of 20" for a direct discharge
- Filters for aspiration and cuspidor bowl are inspected from the outside part

## Spécifications Techniques/Détails :

### Basic unit Partner Art.Nr.SW77020 with Swinging hoses system:

Complete with (Instrumentation From left to right):

- Pos.1=1 3-way f syringe,
- Pos.2=1 Turbine tubing with light
- Pos.3=1 Micro motor tubing with micro motor MC3 LK
- Pos.4=1 Free ( 2<sup>nd</sup> Turbine tubing or micro motor MC3 optional)
- Pos.5=1 Free (EMS or SUPPRASSON Scalar or Curing light Mini LED)
- Control panel for patient's chair and instruments
- Swivel tray support with one tray ( double optional)

### Hydro Group unit :

- Removable and Auto-cleavable ceramic cuspidor bowl
- Can rotate 180° (space saving)
- Cattani cuspidor and Durr selective valves are included
- Safe spray distilled water system

### Assistant table : (Instrumentation from left to right)

- Extensions options to 3 instruments
- Pos.1=1 Surgical section,
- Pos.2=1 saliva ejector,
- Pos.3=Free ( 3-way-syringe for assistant optional)
- The control panel command for the chair functions and instruments

### Fauteuil patient AG04SWP:

- The synchronized trendelenburg movement of the backrest
- 4=3 programs positions, plus additional rising position and return to the last position.
- Left arm rest (Right one optional).
- The Head rest with a double articulation

### The Multifunctional foot control:

- The foot control for the command of all chair's movements and Instruments functions.
- Activation of the instruments with a progressive regulation of micro motors power.
- The spray and chip blower function are included.

### The operating Lamp: (Unit Mounted)

- Equipped with EDI Light (MAIA or ALYA LED Light optional)
- Reliable and proved
- 24.000 Lux, 5000 °k

### Colors :

- For exact definition colors please refer to the separate Chair's upholstery color table
- Unit base, arm systems, and Hydro unit standard color :( color choose optional)

