

AXIOQUICK® ANATOMIC TRANSFERBOW

A light, easy to use anatomic type of facebow referenced to the cranium for more precise positioning of the maxillary cast as it relates to the skull. The interpupillary line leveler, special designed earpieces with hygienic earpiece caps, and nasion relator position the facebow so that a correct patient simulation of the maxillary teeth can be transferred to the articulator. The blue earpiece caps assure a pleasant and comfortable experience for the patient.



TRANSFER FORK ASSEMBLY

Two type of facebow transfer fork assemblies are available. Both use non-torsion type clamps. Transfer fork assemblies can be removed from the anatomic facebow and can be attached to the SAM articulator with facebow transfer stand AX or to most other articulators with other transfer stands.

TRANSFER FORK ASSEMBLY I - NT



- Secure Two-toggle connection
- Removable from facebow
- Non torsion NT clamp mechanism
- Standard facebow transfer fork
- Medical grade steel transfer fork (sterilizable)

TRANSFER FORK ASSEMBLY AX - AXIOMATIC®

- Secure One-toggle non torsion clamp
- Removable from facebow
- · Quick and easy clamping
- Separately removable facebow transfer fork
- Medical grade steel transfer fork (sterilizable)



ATB 350K Transferbow Kit III - NT



ATB 390K Transferbow Kit - AX



Hygienic Earpiece Caps (200 pieces)



ATB 305
Transfer Fork Assembly I-NT



ATB 395
Transfer Fork Assembly AX AXIOMATIC®



ART 398Transfer Stand AX